

Quality Systems Manual

Materials Division

Section 3 -- Equipment Inventory and Calibration

| | |
|---|---------------|
| Inventory and Calibration Information | 3 - 1 |
| Aggregate Lab | 3 - 2 |
| Asphalt Labs | 3 - 5 |
| Bituminous Lab | 3 - 5 |
| Chemistry Lab | 3 - 8 |
| Cement Lab | 3 - 9 |
| Structures and Concrete Lab | 3 - 13 |
| Equipment Calibration Procedures | 3 - 14 |
| Air Meters (cp-am) | 3 - 14 |
| Brass Rings and Assembly (cp-br) | 3 - 16 |
| CBR Annular and Slotted Weights (cp-cbr) and Penetration Piston | 3 - 17 |
| Compression or Loading Device (cp-ct) | 3 - 18 |
| Compressive Strength of Cylindrical Concrete | |
| Specimens Using Neoprene Caps | 3 - 19 |
| Compressive Testing Machine (cp-ct) | 3 - 20 |
| Concrete Test Molds (plastic) (cp-m) | 3 - 21 |
| Conicle Molds and Tamper (cp-cm) | 3 - 22 |
| Direct Tension Machine (cp-dtt) AASHTO T-314 | 3 - 23 |
| Ductility Apparatus (cp-da) | 3 - 26 |
| Dynamic Shear Rheometer (DSR) | 3 - 27 |
| Flash Cleveland Open-Cup (cp-fc) | 3 - 28 |
| Flash Point Tag Open-Cup (cp-fp) | 3 - 29 |
| Flow Gauges (cp-fg) | 3 - 30 |
| Kinematic and Absolute Viscosity Bath (cp- kvis) | 3 - 31 |
| LA Abrasion Machine (cp-la) | 3 - 32 |
| Liquid Limit Device and Grooving Tool (cp-llg) | 3 - 33 |
| Manual Compaction Hammer (cp-mh) | 3 - 34 |
| Marshall Hammer (cp-marh) | 3 - 35 |
| Marshall Mechanical Compactor (cp-mmec) | 3 - 36 |
| Marshall Stability & Flow Test Apparatus (cp-mt) | 3 - 37 |

| | |
|---|-------------------|
| Mechanical Compactor (cp-mc) | 3 - 39 |
| Mechanical Shaker (cp-ms) | 3 - 40 |
| Moist Cabinets, Rooms | 3 - 41 |
| Molds - Soil Testing | 3 - 42 |
| Oven, General Purpose (cp-gpoven) | 3 - 43 |
| Oven, Rolling Thin Film (cp-rtfo) | 3 - 44 |
| Penetrometer (cp-pa) | 3 - 45 |
| Pycnometers (cp-py) | 3 - 46 |
| Rice Flasks (cp-rf) | 3 - 47 |
| Scales and Balances (cp-sb) | 3 - 48 |
| Sieves | 3 - 49 |
| Sieve - 325 | 3 - 50 |
| Slump Cones (cp-sc) | 3 - 52 |
| Sodium Sulfate Containers, Sulfate Oven (cp-ssc) | 3 - 53 |
| Soil Hydrometers (cp-sh) | 3 - 54 |
| Straight Edge | 3 - 55 |
| Thermometers, Liquid-in-Glass (cp-th) | 3 - 56 |
| Timers (cp-ti) | 3 - 57 |
| Unit Weight Container - Concrete (cp-uwct) | 3 - 58 |
| Unit Weight Measures - Aggregate (cp-uwa) | 3 - 59 |
| Vacuum System (cp-vac) | 3 - 60 |
| Water Bath - Penetrometer (cp-wb) | 3 - 62 |
| Calibration Record Forms | 3 - 64 |
| Air Meter (cp-am) | 3 - 64 |
| Brass Rings and Assembly (cp-br) | 3 - 66 |
| CBR (cp-cbr) | 3 - 67 |
| Compression or Loading Device (cp-ct) | 3 - 68 |
| Compressive Test Machine (cp-ct) | 3 - 69 |
| Concrete Test Mold (cp-mcc) | 3 - 71 |
| Conical Molds and Tamper (cp-cm) | 3 - 72 |
| Direct Tension Machine Model | 3 - 73 |
| Dynamic Shear Rheometer (SHRP) (cp-dsr) | 3 - 74 |
| Ductility Apparatus (cp-da) | 3 - 77 |
| Flash Cleveland Open-Cup (cp-fc) | 3 - 79 |
| Flash Point Tag Open-Cup (cp-fp) | 3 - 80 |
| Flow Gauges (cp-fg) | 3 - 81 |
| Kinematic Viscometers | 3 - 82 |
| Liquid Limit Device and Grooving Tool (cp-llg) | 3 - 85 |
| Manual Compaction Hammer (cp-mh) | 3 - 87 |
| Marshall Hammer (cp-marh) | 3 - 88 |
| Marshall Molds & Breaking Head (cp-mmbh) | 3 - 90 |
| Mechanical Shaker (cp-ms) | 3 - 93 |
| Moist Cabinet, Room (cp-mcr) | 3 - 94 |

| | |
|---|----------------|
| Mold -soils testing, compaction (cp-m) | 3 - 95 |
| Oven, General Purpose (cp-gpoven) | 3 - 96 |
| Oven, Rolling Thin Film (cp-rtfo) | 3 - 97 |
| Penetrometer (cp-pa) | 3 - 98 |
| Pycnometers (cp-py) | 3 - 99 |
| Rice Flasks (cp-rf) | 3 - 100 |
| Sand Equivalent Equipment | 3 - 101 |
| Saybolt Viscometer | 3 - 102 |
| Scales and Balances (cp-sb) | 3 - 103 |
| Sieve 325 (cp-sv325)) | 3 - 105 |
| Slump Cone (cp-sc) | 3 - 106 |
| Sodium Sulfate Containers (cp-ssc), | 3 - 107 |
| Soil Hydrometer (cp-sh) | 3 - 108 |
| Straight Edge (cp-se) | 3 - 109 |
| Thermometers (cp-th) | 3 - 110 |
| Thermometer, Liquid-in-Glass (cp-th) | 3 - 112 |
| Timers (cp-ti) | 3 - 113 |
| Unit Weight Container - Concrete (cp-uwct) | 3 - 114 |
| Unit Weight Measures - Aggregate (cp-uwa) | 3 - 115 |
| Vacuum Capillary Viscometers | 3 - 116 |
| Water Bath (cp-wb) | 3 - 119 |
| Water Bath - Penetrometer (cp-wb) | 3 - 120 |
| Weighing Devices - Materials Tests | 3 - 121 |

Inventory and Calibration Information

All Labs

Test Records

- < Each lab supervisor ensures that their equipment is calibrated at the specified intervals
- < calibration record forms are kept in a filing system accessible to the Internal Quality System Reviewer.
- < Any discrepancies in the system are referred to the Engineer for Materials for review and resolution.

Testing

- < All lab test equipment is calibrated on a regular basis.
- < The test records are maintained by each lab supervisor and stored in each lab office.
- < New equipment is calibrated before first use.
- < Equipment that is out of calibration is either repaired, or removed from service.
- < Repaired equipment is calibrated before it is returned to service.

Tracking

Manual card filing system

- < Information is stored in a card filing box with monthly divisions containing one card per piece of equipment. Each card indicates identification of the equipment and its calibration interval.
- < The cards are placed in the month corresponding to the upcoming equipment calibration due date. After each time the equipment is calibrated, the corresponding equipment card is dated and placed in the monthly division of its next calibration due date .
- < The equipment card file is manually checked on a regular basis to ensure all tests are done on time.

Automatic electronic filing system

- < Information is stored electronically using a computer program that combines database functions with automated calendar tracking functions.
- < The computer program stores data about each piece of equipment, its calibration, and its appropriate calibration interval.
- < The computer program notifies the lab supervisor which equipment is due for calibration prior to the calibration date, and automatically resets for the next calibration date.

The spreadsheets on the following pages contain equipment and calibration information for each laboratory.

Aggregate Lab

Accreditation and Certification Information Aggregate Laboratory

| Equipment Item | Manufacturer | Model # | Serial # | UDOT # | Cal./Verif. Proc.Ref | Inter Month | Last Cal. Calibration | Next Calibration | Calibration Person | Date Rec'd | Date in Service | Condition when Rec'd |
|--|------------------|------------|-------------|----------------|----------------------|-------------|-----------------------|------------------|--------------------|------------|-----------------|----------------------|
| Lab Scale 1 | Mettler - Toledo | SG 16001 | 21155097383 | | UDOT cp-sb | 12 | 08/02 | 08/03 | Shafer/Raleigh | 1996 | 1996 | New |
| Lab Scale 2 | Mettler | PM 16 | G01273 | 44-0925 | UDOT cp-sb | 12 | 08/02 | 08/03 | Shafer/Raleigh | 1987 | 1987 | New |
| Lab Scale 3 | Mettler - Toledo | SB 32001 | 1119230904 | | UDOT cp-sb | 12 | 08/02 | 08/03 | Shafer/Raleigh | 2000 | 2000 | New |
| Lab Oven 1 | VWR | 1675 | 900299 | | UDOT cp-do | 4 | 02/03 | 04/03 | Robert Tripp | 1999 | 1999 | New |
| Lab Oven 2 | VWR | 1675 | 1100202 | | UDOT cp-do | 4 | 02/03 | 04/03 | Robert Tripp | 2003 | 2003 | New |
| Lab Oven 3 | VWR | 1685 | 800102 | 44-1369 | UDOT cp-do | 4 | 02/03 | 04/03 | Robert Tripp | 2003 | 2003 | New |
| Mech. Shaker | Rainhart | Mary Ann | | | UDOT cp-ms | 12 | 09/02 | 09/03 | Robert Tripp | 2002 | 2002 | New |
| Mech. Sample Splitter | Gilson | SP - 1 | | 59-40-602 | | | | | | 1970 | 1970 | New |
| Mechanical Sand Equivalent | Soil Test | CL-232 | 144 | 59-44-087 | AASHTO T-17 | 12 | 06/02 | 06/03 | Robert Tripp | 1970 | 1970 | New |
| Sand Equivalent Test Set | Soil Test | CL-230 | | | AASHTO T-17 | 12 | 06/02 | 06/03 | Robert Tripp | 1970 | 1970 | New |
| Fine Void Content Apparatus Test Set | Gilson | SG - 40 | | | | 12 | | | Robert Tripp | 2000 | 2000 | Used |
| Coarse Void Content Apparatus Test Set | Gilson | SG - 42 | | | | 12 | N/A | N/A | Not in Service | 2001 | 2001 | Used |
| Sulfate Soundness Soak Tank | U.D.O.T. | Home made | | | | | N/A | N/A | Not in Service | 1960 | 1960 | New |
| Organic Impurities Test Set | Humboldt | HM-137 | | | | | | | | 1990 | 1990 | New |
| Conical Mold and Tamper | Humboldt | H-3360 | | C - 1 C - 2 | | 12 | 07/02 | 07/03 | Robert Tripp | 1998 | 1998 | New |
| Specific Gravity Bench and Tank | Gilson | SG - 20 | | | | | | | | 2003 | 2003 | New |
| Los Angeles Abrasion Mach. | Soil Test | H-501 | | 59-44-78 | | 24 | 05/01 | 05/03 | Robert Tripp | 1958 | 1958 | New |
| Weight Measures 0.10 | | | | Agg.1 | UDOT cp-uwe | 12 | 05/02 | 05/03 | Robert Tripp | 1997 | 1997 | New |
| 0.10 | | | | Agg.2 | UDOT cp-uwe | 12 | 05/02 | 05/03 | Robert Tripp | 2002 | 2002 | Used |
| 0.25 | Humboldt | H-3664.1 | | Agg.5 | UDOT cp-uwe | 12 | 05/02 | 05/03 | Robert Tripp | 2002 | 2002 | Used |
| 0.50 | Humboldt | H-3661 | | Agg.3 | UDOT cp-uwe | 12 | 05/02 | 05/03 | Robert Tripp | 2002 | 2002 | Used |
| 0.50 | Humboldt | H-3661 | | Agg.4 | UDOT cp-uwe | 12 | 05/02 | 05/03 | Robert Tripp | 2002 | 2002 | Used |
| UA Pulverizer | VWR | 24253-000 | | 44-0945 | | | | | | 1988 | 1988 | New |
| Thermometer | VWR | 8C7992 | | | | 12 | 12/02 | 12/03 | Carter | 1998 | 1998 | New |
| Thermometer | ErTco | L 95 - 781 | | | | 12 | 12/02 | 12/03 | Carter | 1998 | 1998 | New |
| Thermometer | ErTco | 5042 | | | | 12 | 12/02 | 12/03 | Carter | 2001 | 2001 | New |
| Thermometer | ErTco | G 96 - 182 | | | | 12 | 12/02 | 12/03 | Carter | 2000 | 2000 | New |
| Thermometer | | | | | | 12 | 12/02 | 12/03 | Carter | 2000 | 2000 | New |
| Water Bath | Humboldt | | | | | N/A | N/A | N/A | | 2002 | 2002 | Used |
| Sulfate Soundness Containers | UDOT | Home made | | | | | | | | 1997 | 1997 | Used |

| Equipment Item | Manufacturer | UDOT # | Cal./Verif. Proc. Ref | Inter. Month | Last Cal. Calibration | Next Calibration | Calibration Person | Date Rec'd | Date in Service | Condition when Rec'd |
|----------------|--------------|-----------|-----------------------|--------------|-----------------------|------------------|--------------------|------------|-----------------|----------------------|
| 12" SIEVES | | | | | | | | | | |
| 2 1/2 | Endecotts | X000645 | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 1 1/4 | Endecotts | X100875 | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 1 | Tyler | X100945 | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 1 | Endecotts | X100861 | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 1 | Endecotts | X100860 | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 3/4 | Tyler | X001616 | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 3/4 | Endecotts | X100862 | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 3/4 | Endecotts | X100859 | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 5/8 | Tyler | X000884 | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 1/2 | Tyler | X000707A | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 1/2 | Endecotts | X100864 | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 1/2 | Endecotts | X100857 | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 3/8 | Tyler | X000708A | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 3/8 | Endecotts | X100864 | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 3/8 | Endecotts | X100857 | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 5/16 | Fisher | X001003Z | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 1/4 | Tyler | X000690A | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 4 | Tyler | X000601AG | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 4 | Endecotts | X100865 | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 4 | Endecotts | X100856 | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 5 | Endecotts | X100949 | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 8 | Dual | X000638 | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 8 | Tyler | X002638A | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 8 | Endecotts | X100866 | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 8 | Endecotts | X100855 | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 12 | Endecotts | X002926Z | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 12 | Tyler | X000709A | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 16 | Dual | X100180 | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 16 | Endecotts | X100854 | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 16 | Endecotts | X100867 | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 30 | Endecotts | X100868 | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 30 | Endecotts | X100853 | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 40 | Tyler | X0001660B | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 50 | Tyler | X002794Z | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 50 | Endecotts | X100852 | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 50 | Endecotts | X100869 | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 100 | Tyler | X0021137 | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |
| 100 | Endecotts | X100870 | UDOT cp-sd | 6 | 12/02 | 06/03 | Shafer/Raleigh | 1997 | 1997 | New |

| Equipment Item | Manufacturer | UDOT # | Cal./Verif. Proc.Ref | Inter Month | Last Cal. Calibration | Next Calibration | Calibration Person | Date Rec'd | Date in Service | Condition when Rec'd |
|----------------|---------------|----------|-------------------------|----------------|--------------------------|---------------------|-----------------------|---------------|--------------------|----------------------------|
| 12" SIEVES | | | | | | | | | | |
| 100 | Endecotts | X100851 | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 200 | Endecotts | X100874 | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 200 | Tyler | X002139A | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 200 | Endecotts | X101051 | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 200 | Endecotts | X100871 | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 200 | Tyler | X002120 | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 8" SIEVES | | | | | | | | | | |
| 2 1/2 | Tyler | X000678 | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 2 | Tyler | X100923 | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 1 1/2 | Tyler | X00686 | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 1 1/4 | Tyler | X000679 | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 1 | Endecotts | X002431A | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 3/4 | Sargent-Welch | X000648 | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 3/4 | Tyler | X000711A | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 5/8 | Testlab | X100924 | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 5/8 | Testlab | X000676A | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 1/2 | Tyler | X101322 | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 1/2 | Endecotts | X002425 | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 3/8 | Endecotts | X101156 | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 5/16 | Endecotts | X100873 | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 5/16 | Tyler | X000675A | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 1/4 | Tyler | X101321 | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 1/4 | VWR | X100872 | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 4 | Dual | X000629 | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 4 | Tyler | X000604 | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 5 | Tyler | X000674 | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 8 | Gilson | X100199 | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 10 | Tyler | X101282 | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 16 | Fisher | X002169 | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 20 | Tyler | X000713 | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 30 | Endecotts | X101155 | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 40 | Tyler | X000715 | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 50 | Endecotts | X101237 | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 50 | Endecotts | X101236 | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 70 | Tyler | X100926 | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 80 | Tyler | X000626 | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 100 | Tyler | X000717 | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 100 | Tyler | X101159 | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |
| 200 | Dual | X002198 | UDOT cp-sd | 6 | 12/02 | 12/02 | Shafer/Raleigh | 1997 | 1997 | New |

Asphalt Labs

Asphalt Laboratory

Date: 01/14/2003

| Equipment | Date Rec'd | Date in Service | Manufacturer | Model Number | Serial Number | Condition received |
|------------------------------|------------|-----------------|----------------------|------------------|---------------------|--------------------|
| Rolling-Thin Film Oven | 7/83 | 7/83 | Cox and Sons | CS-325-A | 44-717 | new |
| Rolling-Thin Film Oven | 10/93 | 10/93 | Cox and Sons | CS-325-A | 44-1177 | used |
| Rolling-Thin Film Oven | 6/2000 | 2/2001 | Cox and Sons | CS-325-B | 44-1335 | new |
| Laboratory Oven | 12/21/99 | 12/14/99 | VWR | 1350GM | 1100302 | new |
| Laboratory Oven | 10/91 | 10/91 | Lab Line | Imperial V 3486 | 44-1168 | new |
| Laboratory Oven | 2/93 | 2/93 | Lab Line | Imperial V 3495 | 0293-0417 | new |
| Viscometer Bath | 3/95 | 3/95 | Cannon | CT-1000 | 44-1203 | new |
| Viscometer Bath | 2-01 | 2-01 | Cannon | CT-1000 | 44-1164 | used |
| Saybolt Viscometer Bath | 8/92 | 8/02 | Precision Scientific | 74966 | 44-573 | used |
| Saybolt Viscometer Bath | 1/2001 | 6-2001 | Koehler | K21410 | R90190065 | new |
| Vacuum Regulator | 2/95 | 2/95 | Cannon | NA | 44-890 | new |
| Rotational Viscometer | 1/93 | 1/93 | Brookfield | RV-DV-11+ | DE25287 | new |
| Dynamic Shear Rheometer | 1/95 | 1/95 | B ohlin | DSR 11 | 95/3520/BOH/01a/008 | new |
| Dynamic Shear Rheometer | Sept 98 | Sept 98 | Bohlin | CVO 50 | 98006670curc112 | new |
| Bending Beam Rheometer | Feb 2002 | March02 | Cannon | TE-BBR | 3294-A502 | new |
| Bending Beam Rheometer | Apr 98 | Apr 98 | Cannon | TEBBR | 3127-C | new |
| Pressure Aging Mechanism | July 99 | July 99 | Prentex | 9300 | 93154 | new |
| Direct Tension Machine | 9/1999 | 9/1999 | Instron Asphalt Pro | 5525 | 5525-C9827 | new |
| Direct Tension Machine | 1-97 | 1-97 | Instron BTI 3 | 55BT1 | C6646 | new |
| Cleveland Open-Cup | 7/90 | 7/90 | Humboldt | gas | 44-0001 | new |
| Ductility Machine | 1-1994 | 1-1994 | Forma Scientific | 2816 | 20457-16 | new |
| Penetrometer | 7/85 | 7/85 | Precision Scientific | NA | 44-851 | new |
| Water Bath | 1/91 | 1/91 | Forma | 2095 circulating | 20974-092 | new |
| Water Bath | 1/91 | 1/91 | 2095 Forma | 2095 cirulating | 20962-35 | new |
| Water Bath | 1/91 | 1/91 | Forma | circulating | 20154-1392 | new |
| Linear measuring: | | | | | | |
| Caliper | 3/97 | 3/97 | Mitutoyo | CD-6" CS | 50312 | new |
| NIST Traceable Thermometers: | | | | | | |
| ASTM 45F | 94-95 | 94-95 | Brooklyn | liquid in glass | 2N359 | new |
| minus 30 to 2C(62C) | 94-95 | 94-95 | Princo | liquid in glass | BBB09619 | new |
| ASTM 63F | 94-95 | 94-95 | Brooklyn | liquid in glass | 9161399 | new |
| ASTM 96C | 3-2001 | 3-2001 | VWR | liquid in glass | 3L0644 | new |
| ASTM 13C | 94-95 | 94-95 | Brooklyn | liquid in glass | 1D622 | new |
| ASTM 19C | 94-95 | 94-95 | Brooklyn | liquid in glass | 33849 | new |
| ASTM 47F | 94-95 | 94-95 | Brooklyn | liquid in glass | 54931 | new |
| ASTM 110F | 12/96 | 12/96 | Princo | liquid in glass | 689810 | new |
| ASTM 94C | 2-200 | 2-2000 | A&M | liquid in glass | 2024180 | new |
| SERIES II MULTIMETER | 8-2000 | 8-2000 | Fluke | 77 | 62610919 | new |
| TURE RMS MULTIMETER | 6-1999 | 6-1999 | Fluke | 8060A | 663505 | new |
| Pressure/Vaccum Transducer | 6-1999 | 6-1999 | Fluke | PV350 | id# 2151A32 | new |
| Rheometer Temp Probe | 6-2001 | 6-2001 | Cannon | 9728-V95 | DSR 100K/2246 | new |
| Rheometer Temp Probe | 10-1996 | 10-1999 | Cannon | 9728-V95 | DSR 100K/1497 | new |
| Electronic balance | 7/86 | 7/86 | Sartorius | 1465MP8-2 | 44-873 | new |
| Electronic Balance | 9/93 | 9/93 | Mettler | AE200S | 44-1171 | new |
| Electronic Balance | 9/91 | 9/91 | Mettler | PM-11N | 44-1118 | new |

Note: Manufacturer's instructions may be found in this lab supervisor's office.

Bituminous Lab

| Equipment Calibration and Verification Information BITUMINOUS LAB | | | | | | | | | | | | |
|--|-------------------|------------|----------|---------|----------------|-------------------|------------------|------------------|----------------------------|------------|-----------------|-----------------|
| Equipment Item | Manufacturer | Model # | Serial # | UDOT # | Calib/Verif | Interval (Months) | Last Calibration | Next Calibration | Calibration Responsibility | Date Rec'd | Date In-Service | Condition Rec'd |
| | | | | | Procedure Ref. | | | | | | | |
| Analytic Balance | Mettler | PC-16 | D42362 | 44-0637 | UDOT | 12 | Jan-03 | Jan-04 | Starkie | 1984 | 1984 | New |
| Analytic Balance | Mettler | PE-16 | D45600 | 44-864 | UDOT | 12 | Jan-03 | Jan-04 | Starkie | 1988 | 1988 | New |
| Analytic Balance | Mettler | PM-11 | G15383 | 44-0924 | UDOT | 12 | Jan-03 | Jan-04 | Starkie | 1991 | 1991 | New |
| Analytic Balance | Mettler | PM-3000K | 21134658 | | UDOT | 12 | Jan-03 | Jan-04 | Starkie | 1994 | 1994 | New |
| Vacuum System | Welch Duo-Seal | 1402 L-05 | 130542 | 44-1122 | UDOT | 6 | Jan-03 | Aug-03 | Steve N. | 1993 | 1993 | New |
| Mechanical Shaker | Gilson | SS-12R | 1004 | 44-615 | UDOT | 12 | Jan-03 | Jan-04 | Steve N. | 1985 | 1985 | New |
| Oven (Gas) | Hussman Southbend | Marathoner | 1165787 | 44-0972 | UDOT | 4 | Jan-03 | Jun-03 | Steve N. | 1988 | 1988 | New |
| Oven (Electric) | VWR Scientific | 1645-D | N/A | 44-1074 | UDOT | 4 | Jan-03 | Jun-03 | Steve N. | 1991 | 1991 | New |
| Oven (Electric) | VWR Scientific | 1685 | N/A | 44-1143 | UDOT | 4 | Jan-03 | Jun-03 | Steve N. | 1993 | 1993 | New |
| Water Bath #1 | VWR Scientific | 1235-PC | 701593 | N/A | UDOT | 6 | Jan-03 | Aug-03 | Steve N. | 1993 | 1993 | New |
| Water Bath #2 | VWR Scientific | 1235-PC | 401193 | N/A | UDOT | 6 | Jan-03 | Aug-03 | Steve N. | 1993 | 1993 | New |
| Water Bath #3 | VWR Scientific | 1235-PC | 202993 | N/A | UDOT | 6 | Jul-03 | Jan-04 | Steve N. | 1993 | 1993 | New |
| Sieve 3/4" | Tyler | N/A | 001323B | N/A | UDOT | 6 | Oct-02 | Apr-03 | G. Carter | 1990 | 1990 | New |
| Sieve 1/2" | Tyler | N/A | 000134B | N/A | UDOT | 6 | Oct-02 | Apr-03 | G. Carter | 1990 | 1990 | New |
| Sieve 3/8" | Tyler | N/A | 000116B | N/A | UDOT | 6 | Oct-02 | Apr-03 | G. Carter | 1990 | 1990 | New |
| Sieve No. 4 | Tyler | N/A | 002453B | N/A | UDOT | 6 | Oct-02 | Apr-03 | G. Carter | 1990 | 1990 | New |
| Sieve No. 8 | Tyler | N/A | 002457B | N/A | UDOT | 6 | Oct-02 | Apr-03 | G. Carter | 1990 | 1990 | New |
| Sieve No. 16 | Tyler | N/A | 002454B | N/A | UDOT | 6 | Oct-02 | Apr-03 | G. Carter | 1990 | 1990 | New |
| Sieve No. 30 | Tyler | N/A | 000163B | N/A | UDOT | 6 | Oct-02 | Apr-03 | G. Carter | 1990 | 1990 | New |
| Sieve No. 50 | Tyler | N/A | 000042B | N/A | UDOT | 6 | Oct-02 | Apr-03 | G. Carter | 1990 | 1990 | New |
| Sieve No. 100 | Tyler | N/A | 000169B | N/A | UDOT | 6 | Oct-02 | Apr-03 | G. Carter | 1990 | 1990 | New |
| Sieve No. 200 | Tyler | N/A | 002855B | N/A | UDOT | 6 | Oct-02 | Apr-03 | G. Carter | 1990 | 1990 | New |
| Sieve No. 16 (Washing Only) | Tyler | N/A | 002850B | N/A | UDOT | 6 | Oct-02 | Apr-03 | G. Carter | 1990 | 1990 | New |
| Sieve No. 200 (Washing Only) | Tyler | N/A | 000167B | N/A | UDOT | 6 | Oct-02 | Apr-03 | G. Carter | 1990 | 1990 | New |
| Lab Timer #1 | Cole-Parmer | N/A | 21383236 | N/A | UDOT | 6 | Jul-03 | Jan-04 | Steve N. | 2002 | 2002 | New |
| Lab Timer #2 | Cole-Parmer | N/A | 21383244 | N/A | UDOT | 6 | Jul-03 | Jan-04 | Steve N. | 2002 | 2002 | New |
| Thermometer | Ertco | 2F BF | 96-180 | N/A | UDOT | 6 | Jan-03 | Aug-03 | Steve N. | 1994 | 1994 | New |
| Thermometer | Ertco | 2F BF | 96-158 | N/A | UDOT | 6 | Jan-03 | Aug-03 | Steve N. | 1994 | 1994 | New |
| Thermometer | VWR | 2F | 3242 | N/A | UDOT | 6 | Jan-03 | Aug-03 | Steve N. | 1994 | 1994 | New |
| Thermometer | Ertco | 2F | C 96-116 | N/A | UDOT | 6 | Jan-03 | Aug-03 | Steve N. | 1994 | 1994 | New |
| Thermometer | Ertco | 7F | L 95-800 | N/A | UDOT | 6 | Jan-03 | Aug-03 | Steve N. | 1994 | 1994 | New |
| Thermometer | Fisher | 7F | IV9396 | N/A | UDOT | 6 | Jan-03 | Aug-03 | Steve N. | 1994 | 1994 | New |
| Thermometer | Brooklyn | 7F | 69710 | N/A | UDOT | 6 | Jan-03 | Aug-03 | Steve N. | 1994 | 1994 | New |
| Thermometer | Fisher | 7F | 2B5793 | N/A | UDOT | 6 | Jan-03 | Aug-03 | Steve N. | 1994 | 1994 | New |
| Thermometer | Ertco | 7F | K 95-641 | N/A | UDOT | 6 | Jan-03 | Aug-03 | Steve N. | 1994 | 1994 | New |
| Thermometer | Fisher | 17F | IS5898 | N/A | UDOT | 6 | Jan-03 | Aug-03 | Steve N. | 1994 | 1994 | New |
| Thermometer | Princo | 17F | A00421 | N/A | UDOT | 6 | Jan-03 | Aug-03 | Steve N. | 1994 | 1994 | New |
| Thermometer | Ertco | 47F | 1137 | N/A | UDOT | 6 | Jan-03 | Aug-03 | Steve N. | 1994 | 1994 | New |

| | | | | | | | | | | | | |
|---|-----------------------|---------------|--------------|---------|------|----|--------|--------|------------|------|------|-----|
| Flow Gauge | Gilmont | 2P | N/A | N/A | UDOT | 12 | Jan-03 | Jan-04 | Steve N. | 1992 | 1992 | New |
| Superpave Gyratory Compactor | Troxler | 4140 | 203 | N/A | UDOT | 6 | Jan-03 | Aug-03 | Steve N. | 1995 | 1995 | New |
| Superpave Gyratory Compactor | Test Quip | BSCG 1 | | N/A | UDOT | 6 | Jan-03 | Aug-03 | Steve N. | 1995 | 1995 | New |
| Oven Electric | VWR SCIENTI FIC | 1685 | 700398 | 44-1305 | UDOT | 4 | Jan-03 | Jun-03 | Steve N. | 1999 | 1999 | New |
| Ignition Oven 1 | BS Thermal Line | | | | UDOT | 12 | Jan-03 | Jan-04 | G. Starkie | 1991 | 1991 | New |
| Ignition Oven 2 | BS Thermal Line | | | | UDOT | 12 | Jan-03 | Jan-04 | G. Starkie | 1991 | 1991 | New |
| Eletronic Thermometer | Fluke (on APA) | Model 52 | 6452269 | N/A | UDOT | 6 | Jan-03 | Aug-03 | G.Carter | 1994 | 1994 | New |
| Eletronic Thermometer | Fluke(Ov en) | Model 52 | AA00053 7 | N/A | UDOT | 6 | Jan-03 | Aug-03 | G. Carter | 1994 | 1994 | New |
| Eletronic Thermometer | Fluke(Ov en) | Model 52II | 7610028 7 | N/A | UDOT | 6 | Jan-03 | Aug-03 | G.Carter | 2002 | 2002 | New |
| Superpave Gyratory Compactor | Pine | AFGB1 A | 5195 | | UDOT | 6 | Feb-03 | Sep-03 | Steve N. | | | |
| Note: Calibration and verification records may be found in this lab supervisor's office. | | | | | | | | | | | | |

Chemistry Lab Inventory

Utah Department of Transportation: Materials Division Inventory of Major Equipment and Reference Materials

Chemistry Lab

Date: 12/04/2002

| Equipment | Date Rec'd | Date in Service | Manufacturer | Model Number | Serial Number | UDOT Number | Condition when received |
|---------------------------------|----------------------------|-----------------|-----------------------|----------------|----------------|-------------|-------------------------|
| Infrared Spectrophotometer | 1999 | 07/09/99 | Shimadzu | 8300 | A20913600758LP | 44-1312 | new |
| Mass Spectrometer | 0.0024975 02497502 5 | 05/16/02 | Agilent | 5973N | US10442317 | 44-1357 | |
| Gas Chromatograph | 12/30/99 | 05/16/02 | Agilent | 6890N | US10208027 | 44-1357 | new |
| Atomic Absorption Spectrometer | 1991 | 1991 | Perkin Elmer | 3100 | 148202 | 44-1150 | new |
| X-Ray Fluorescence Spectrometer | 5-20-97 | 6-1-97 | Oxford Instruments | MDX-1000 | MX 1142 | 44-1259 | new |
| CertiPrep Mixer/Mill | 5-20-97 | 6-1-97 | SPEX | 8000 | 97023 | 44-1259-01 | new |
| Hydraulic Press | 5-20-97 | 6-1-97 | Carver | C | 41000-341 | 44-1259-02 | new |
| Muffle Furnace - Large | 1980 | not in service | Thermolyne | F-1930-1 | 17201 | 59-40-677 | new |
| Muffle Furnace - Large | 04/11/00 | 04/11/00 | Thermolyne | F48020 | 1058990220624 | | new |
| Muffle Furnace - Large | 04/11/00 | 04/11/00 | Thermolyne | F6000 | 1060000108824 | | new |
| Muffle Furnace - Small | 1980 | 1980 | Thermolyne | FB1415M | 34802762 | 44-869 | new |
| Muffle Furnace - Small | 10-02-96 | 10-02-96 | Thermolyne | FB1415M | 746960706100 | | new |
| Hot Plate | 1980 | 1980 | Thermolyne | HP36025/VWR300 | 36000737 | | new |
| Hot Plate | 1980 | 1980 | Thermolyne | HPA1915B | 06417808 | | new |
| Hot Plate | 1980 | 1980 | Thermolyne | HPA1915B | | | new |
| Hot Plate | 1980 | 1980 | Thermolyne | HPA1915B | | | new |
| Hot/Stir Plate | 1980 | 1980 | Thermolyne | SP - A1025B | 670165 | | new |
| Hot/Stir Plate | 04/23/01 | 04/23/01 | VWR | 357 | 2054 | | new |
| Incubator | 12/05/98 | 12/08/98 | Yamato | IC600 | A6700039 | | new |
| Drying Oven | 1980 | 1980 | Lab-Line Instruments | | 0666 | 40-505 | new |
| Drying Oven | 1980 | 1980 | Blue M Electric Co. | OV-185A | RP 416 | 59-40-645 | new |
| Water/Oil Bath | 1975 | not in service | Lab-Line Instruments | 3005-7 | 0269 | 59-40-288 | new |
| Water/Oil Bath | 1975 | 1975 | Lab-Line Instruments | 3005-7 | 1065 | | new |
| pH Meter | 1975 | 1975 | Sargent-Welch | S-29999/NX | 2509070 | 59-40-384 | new |
| Analytical Balance | 1993 | 1993 | Mettler | AE163 | C-01822 | 44-707 | new |
| Top-loader Balance | 03/18/02 | 03/18/02 | Mettler Toledo | PB602S | 1121082579 | | new |
| Pulverizer | 5-88 | 5-88 | BICO | UA53 | 67463 | 44-0945 | new |
| Vacuum Pump | 3-24-98 | 3-24-98 | Gast | SA55NXGTE-4870 | 9710714031 | | new |
| Refridgerator | 1980 | 1980 | Frigitemp Corp | SS-65-ASCRTM | 4110010078153 | 59-40-27 | new |
| TCLP Extractor | 1993 | 1993 | Anal. Testing & Cons. | | 8054418 | 44-1132 | new |
| Hot Plate | 3-16-98 | 3-16-98 | Thermolyne/Cimarec | HP47135 | 1073971142582 | | new |
| Hot Plate | 3-16-98 | 3-16-98 | Thermolyne/Cimarec | HP47135 | 1073971142575 | | new |
| Spectrophotometer | 7/98 | 7/98 | Hach | DR20/10 | 980500008868 | 44-1291 | new |
| Ion Chromatography system | 06/08/98 | 07/12/98 | Lachat/Zellweger | QuickChem 8000 | A83000-1224 | 44-1290 | new |
| Reagent Pump | 06/08/98 | 07/12/98 | Lachat/Zellweger | IC-110V | A8200-495 | 44-1290-02 | new |
| Eluent Pump | 06/08/98 | 07/12/98 | Lachat/Zellweger | 241485 | A28911-455 | 44-1290-01 | new |
| Top-loader Balance | pre-1978 | not in service | Mettler | K7T | G-1236 | 59-40-562 | new |
| Top-loader Balance | pre-1989 | not in service | Mettler | P160N | 313783 | 59-40-274 | new |
| Analytical Balance | pre-1989 | not in service | Ainsworth | Type 10 | 43635 | 59-40-27 | new |

Note: Manufacturer's instructions may be found in this lab supervisor's office.

Cement Lab Inventory

| Equipment Item | Manufacture | Model # | Serial # | UDOT # | Calibration Procedure Ref. | Interval (months) | Last Calibration | Next Calibration | Calibration Responsibility |
|-----------------|---------------|----------|------------|---------|----------------------------|-------------------|------------------|------------------|----------------------------|
| QSM | UDOT | | | | | | | | T. Peterson |
| Balances: | | | | | | | | | |
| Analytical | Mettler | PM 1200 | SRN-J66801 | 44-1024 | UDOT cp-sb | 12 | 8/02 | 8/03 | G. Starkie |
| Gen. Purpose | Mettler | PC 16 | D42362 | 44-0637 | UDOT cp-sb | 12 | 12/02 | 12/03 | G. Starkie |
| Thermometers: | | | | | | | | | |
| Glass - 300 F | Sargent | ASTM 1F | 43952 | | UDOT cp-th | 6 | 10/02 | 4/03 | G. Carter |
| Glass - 300 F | VWR | ASTM 1F | 600875 | | UDOT cp-th | 6 | 10/02 | 4/03 | G. Carter |
| Glass - 580 F | ERTCO | ASTM 2F | G77737 | | UDOT cp-th | 6 | 10/02 | 4/03 | G. Carter |
| Glass - 305 C | ERTCO | ASTM 69C | 1901 | | UDOT cp-th | 6 | 10/02 | 4/03 | G. Carter |
| Glass - 80 F | VWR | ASTM 17F | A01604 | | UDOT cp-th | 6 | 10/02 | 4/03 | G. Carter |
| Sieves: | | | | | | | | | |
| #20, 12" Screen | Tyler | | X002016 | | UDOT cp-sv | 6 | 10/02 | 4/03 | G. Carter |
| #20, 12" Screen | Tyler | | X100171 | | UDOT cp-sv | 6 | 10/02 | 4/03 | G. Carter |
| #20, 12" Screen | Tyler | | X000728 | | UDOT cp-sv | 6 | 10/02 | 4/03 | G. Carter |
| #16, 8" Screen | Fisher Scien. | | X002746 | | UDOT cp-sv | 6 | 10/02 | 4/03 | G. Carter |
| #20, 8" Screen | Fisher Scien. | | X002019 | | UDOT cp-sv | 6 | 10/02 | 4/03 | G. Carter |
| #40, 8" Screen | Fisher Scien. | | X000745 | | UDOT cp-sv | 6 | 10/02 | 4/03 | G. Carter |
| #50, 8" Screen | Fisher Scien. | | X002017 | | UDOT cp-sv | 6 | 10/02 | 4/03 | G. Carter |
| #100, 8" Screen | Fisher Scien. | | X002018 | | UDOT cp-sv | 6 | 10/02 | 4/03 | G. Carter |
| #20, 8" Screen | Endecott | | X100173 | | UDOT cp-sv | 6 | 10/02 | 4/03 | G. Carter |
| #20, 8" Screen | Endecott | | X100172 | | UDOT cp-sv | 6 | 10/02 | 4/03 | G. Carter |
| #30, 8" Screen | Sargent | | X100782 | | UDOT cp-sv | 6 | 10/02 | 4/03 | G. Carter |
| #4, 8" Screen | Sargent | | X100725 | | UDOT cp-sv | 6 | 10/02 | 4/03 | G. Carter |

| | | | | | | | | | |
|----------------------|---------------|-----------|---------|---------|----------------------|--------------|-------|------|-----------|
| #8, 8" Screen | Fisher Scien. | | X000730 | | UDOT cp-sv | 6 | 10/02 | 4/03 | G. Carter |
| 325 sieve cups (1-9) | Humboldt | | | | See ASTM C-430 | After 5 uses | 10/02 | | R. Davis |
| Drying Oven | Blue M | SW-11TA-1 | SW-4376 | | | | | | |
| Autoclave | Boekel | 95-50-11 | | 44-1222 | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| Molds: | | | | | | | | | |
| 2x2" Cube | Humboldt | | Y-3608 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| 2x2" Cube | Humboldt | | Y-3609 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| 2x2" Cube | Humboldt | | Y-3610 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| 2x2" Cube | Humboldt | | Y-3611 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| 2x2" Cube | Humboldt | | Y-3612 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| 2x2" Cube | Humboldt | | Y-3613 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| 2x2" Cube | Humboldt | | Y-3614 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| 2x2" Cube | Humboldt | | Y-3615 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| 2x2" Cube | Humboldt | | Y-3616 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| 2x2" Cube | Humboldt | | Y-3617 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| 2x2" Cube | Humboldt | | Y-3618 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| 2x2" Cube | Humboldt | | Y-3619 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| 2x2" Cube | Humboldt | | Y-3620 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| 2x2" Cube | Humboldt | | Y-3621 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| 2x2" Cube | Humboldt | | Y-3622 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| 2x2" Cube | Humboldt | | Y-3623 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| Cone Molds | Humboldt | | E-7533 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| Cone Molds | Humboldt | | R-749 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| Cone Molds | Humboldt | | R-750 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| Cone Molds | Humboldt | | S-4119 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |

| | | | | | | | | | |
|---------------------|---------------------|----------|----------------|---------|-------------------------|----|-------|-------|------------|
| Cone Molds | Humboldt | | W-2251 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| Cone Molds | Humboldt | | W-2252 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| Cone Molds | Humboldt | | W-2253 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| Cone Molds | Humboldt | | W-2254 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| Cone Molds | Humboldt | | W-2256 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| Cone Molds | Humboldt | | W-2257 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| Cone Molds | Humboldt | | W-2258 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| Cone Molds | Humboldt | | W-2259 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| Cone Molds | Humboldt | | Y-2161 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| Cone Molds | Humboldt | | Y-2162 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| Cone Molds | Humboldt | | Y-2164 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| Cone Molds | Humboldt | | Y-870 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| Autoclave Molds | Soil Test | | | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| Autoclave Molds | Humboldt | 9252 | | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| Moist Cabinet | Forma Scientific | 3911 | 99343-97 | 441326 | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| Blaine Apparatus | Humboldt | | S4200 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| Flow Table | Humboldt | | | 44-0292 | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| Vicat Apparatus | Humboldt | | W2248 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| Mech. Mixer | Hobart | N-50 | 31-1167-616 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| 400 ml. Cup | | | | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| Comp. Machine | Satec | Mark III | 2500GP187 8 | 441330 | UDOT cp-ct | 12 | 12/02 | 12/03 | G. Starkie |
| N.C. Apparatus | Fisher Scien. | | A602 | | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| Comparator | Humboldt | | | 44-0398 | CCRL inspection used | 30 | 12/02 | 6/05 | CCRL |
| No. 325 Nozzle | | | | | See ASTM C- 430 | 6 | 12/02 | 6/03 | R. Davis |

| | | | | | | | | | |
|-----------------|---------|-------------|--------|--|----------------|---------------|-------|------|------------|
| Blaine Standard | NIST | 114p | | | See ASTM C-204 | 6 | 12/02 | 6/03 | R. Davis |
| Bearing Blocks | | | | | UDOT cp-ct | 12 | | | G. Starkie |
| Std. Sand | AcuSand | C-109+20/30 | | | See ASTM C-778 | upon recieval | | | R. Davis |
| Record. Therm. | VWR | ASTM 17F | A01604 | | See ASTM C-511 | 6 | 12/02 | 6/03 | R. Davis |

Structures and Concrete Lab Inventory

| Equipment Item | Manufacturer | Model # | Serial # | UDOT # | Cal./Verif. Proc.Ref | Inter Month | Last Cal. Calibration | Next Calibration | Calibration Person | Date Rec'd | Date in Service | Condition when Rec'd |
|--------------------------------|----------------------|----------------|------------------|---------|-------------------------|----------------|--------------------------|---------------------|-----------------------|---------------|--------------------|-------------------------|
| Large Press | Reihle | FH-450 | R-94251 | | UDOT cp-ct | 12 | 12/02 | 12/03 | Graham Starkie | 1960 | 1960 | New |
| Reinforcing Steel Bender | UDOT | N/A | None | | UDOT | | | | Robert Winters | 1985 | 1985 | New |
| Skidmore Machine | Skidmore- Wilheim | ML | 9136 | | UDOT | | | | Robert Winters | 01/01/81 | 01/01/81 | New |
| Rockwell Hardness Tester | Wilson Instrument | 4YR a | 367 | | UDOT | | | | Robert Winters | 06/01/75 | 06/01/75 | New |
| Large Scale | Neweigh | S20T | 45307 | 40-0181 | UDOT cp-sb | 12 | | | Graham Starkie | 01/01/81 | 01/01/81 | New |
| Extensometer | Tinius Olson | S-8187 | 180976 | | UDOT cp-ct | | | | Robert Winters | | | New |
| Instron Press | Instron | 1000 | 2703 | | UDOT cp-ct | 12 | Out of Service | | Graham Starkie | 03/15/83 | 03/20/85 | New |
| Mullen Tester | BF Perkins | N/A | 8205-84- 1994 | | UDOT | | | | Robert Winters | 1968 | 1968 | New |
| Satec Press | Satec | Apex T10000 | 1064 | | UDOT cp-ct | 12 | Out of Service | | Graham Starkie | 1996 | 1996 | New |
| Shore Durometer | Shore Instrument | A-2 | 70810 | | UDOT | | | | Robert Winters | 1978 | 1978 | New |
| Slump Cone | Humboldt | | | CL-1 | UDOT cp-sc | 12 | 01/03 | 01/04 | Bill Lawrence | 1999 | 1999 | New |
| Slump Cone | Humboldt | | | CL-2 | UDOT cp-sc | 12 | 01/03 | 01/04 | Bill Lawrence | 1999 | 1999 | New |
| Slump Cone | Gilson | HM-45 | | CL-5 | UDOT cp-sc | 12 | 01/03 | 01/04 | Bill Lawrence | 2003 | 2003 | New |
| Concrete Pressure Meter | Watts | | | BZ-05 | UDOT cp- am | 12 | 03/02 | 03/03 | Gerald Carter | 2002 | 2002 | New |
| Concrete Pressure Meter | Watts | | | CL-34 | UDOT cp- am | 12 | 09/02 | 09/03 | Gerald Carter | 2002 | 2002 | New |
| Concrete Pressure Meter | Watts | | | CL-35 | UDOT cp- am | 12 | 03/03 | 03/04 | Gerald Carter | 2002 | 2002 | New |
| Concrete Pressure Meter | Watts | | | CL-60 | UDOT cp- am | 12 | 02/02 | 02/03 | Gerald Carter | 2002 | 2002 | New |
| Concrete Pressure Meter | Watts | | | CL-61 | UDOT cp- am | 12 | 02/02 | 02/03 | Gerald Carter | 2002 | 2002 | New |
| Concrete Pressure Meter | Watts | | | CL-62 | UDOT cp- am | 12 | 12/02 | 12/03 | Gerald Carter | 2002 | 2002 | New |
| Concrete Pressure Meter | Watts | | | CL-63 | UDOT cp- am | 12 | 02/03 | 02/04 | Gerald Carter | 2002 | 2002 | New |
| Concrete Pressure Meter | Watts | | | CL-72 | UDOT cp- am | 12 | 02/03 | 02/04 | Gerald Carter | 2002 | 2002 | New |
| Concrete Pressure Meter | Watts | | | 40-812 | UDOT cp- am | 12 | 02/03 | 02/04 | Gerald Carter | 2002 | 2002 | New |
| Concrete Pressure Meter | Watts | | | 40-821 | UDOT cp- am | 12 | 02/02 | 02/03 | Gerald Carter | 2002 | 2002 | New |

Equipment Calibration Procedures

Air Meters (cp-am)

Type of Equipment: Pressure Meter

AASHTO T 152, C 231

Purpose: To verify the accuracy of a pressure type air meter.

Inspection Equipment

1. Threaded straight tubing
2. Threaded bent tubing
3. Syringe

Procedure

1. Fill the base with water.
2. Screw the short piece of straight tubing into the threaded petcock hole on the underside of the cover.
Clamp cover on the base with the tube extending down into the water.
3. With both petcocks open, add water with a syringe through the petcock having the pipe extension below, until all air is forced out through the opposite petcock. Leave both petcocks open.
4. Pump up air pressure to a little beyond the predetermined initial pressure line. Wait a few seconds for compressed air to cool to normal temperature and then stabilize the gauge hand at the proper initial pressure line by pumping or bleeding off air as needed.
5. Close both petcocks and immediately press down on the thumb lever exhausting air into the base.
Wait a few seconds until the hand is stabilized.
If all the air was eliminated and the initial pressure line was correctly selected, the gauge should read 0%.
If two or more tests show a consistent variation from 0% in the result, then change initial pressure line to compensate for the variation. Use the newly established "initial pressure" line for subsequent tests.
6. Screw curved tube into the outer end of petcock.
Pressing on thumb lever and controlling flow with petcock lever, fill the 5% calibrating vessel (345 ml) level full of water from the base.
7. Release the air at the free petcock. Open the other petcock and let the water in the curved pipe run back into the base. There should be 5% air in the base.
8. With petcocks open, pump air pressure in exact manner as outlined in paragraph 4. Close

petcocks and immediately press the thumb lever. Wait a few seconds for exhaust air to warm to normal temperature, and for the needle to stabilize. The dial should now read 5%.

9. If two or more consistent tests show that the gauge reads more or less than 5% air by .2% (or whatever is considered satisfactory), then remove gauge glass and reset the dial hand to 5% by turning the recalibrating screw located just below and to the right of the center diameter.
10. When gauge hand reads correctly at 5%, additional water may be withdrawn in same manner to check results at 10%, 15%, 20%, and so on.

Brass Rings and Assembly (cp-br)

T-53

Purpose To measure the physical dimensions of the brass ring assembly.

Inspection Equipment Required

1. Balance capable of weighing 2000 g., readable to 0.001 g.
2. Calipers readable to 0.01 mm.

Tolerance Tolerances can be found in the test method listed above.

Procedure

1. With calipers, measure and record the dimensions for the rings, as shown in T-53 fig. 1(a).
2. With calipers, measure and record the dimensions for the ring holder, as shown in T-53 fig.1(b)
3. With calipers, measure and record the dimensions for the ball centering guide as shown in T-53 fig.1 (c).
4. With calipers, measure and record the diameter of the two steel balls.
5. On a balance or scale, weigh each ball and record the weights.

CBR Annular and Slotted Weights (cp-cbr) and Penetration Piston

T193/D1883

Purpose To measure critical dimensions and weights of the CBR Annular and Slotted Weights.

Inspection Equipment Required

1. Ruler readable to the nearest 1/16 inch.
2. Balance, capacity 15 lbs., readable to 0.01 lbs.

Tolerance: CBR annular and slotted weights and penetration piston must meet the dimensional tolerances and be within the specified weight tolerances as specified in the test method.

Procedure

1. Measure and record the dimensions of each weight.
2. Weigh and record the weight of each weight.
3. Measure and record the penetration piston shaft diameter and the penetration piston length.

Compression or Loading Device (cp-ct)

T193/D1883, T208/D2166, T216/D2435, T234/D2850, T236/D3080

Purpose

To check the critical components of compression loading devices used in the above referenced testing methods.

Inspection Equipment Required

A load cell which is calibrated to the accuracy required for the specific compression device and test method.

Tolerance

The equipment will meet the tolerances specified in the test methods referenced above.

Procedure

1. Insert calibrated load cell and compression device where the pressure is applied.
2. Observe and record compression device load indications throughout the range of the loads required for the test methods referenced above according to AASHTO T67.
3. Verify that applications of loads can be applied within the time limits specified in the above referenced test procedures.

Compressive Strength of Cylindrical Concrete Specimens Using Neoprene Caps

T22-92, D2000

Purpose To test compressive strength of cylindrical concrete specimens using neoprene caps instead of sulphur caps.

Inspection Equipment Required

As required in ASTM-D2000 and AASHTO-T 22-92.

Procedure

The procedure followed is detailed in the main test method, (AASHTO T22-92, Sec. 7 and ASTM C39-86), modified as noted below:

1. Place an extrusion controller which contains a neoprene cap on the top and bottom surfaces of the concrete cylinder.
With the neoprene caps in contact with the concrete cylinder, carefully align the axis of the specimen with the center of thrust of the spherically seated block. Bring the bearing blocks of the machine in contact with both of the extrusion controllers.
2. Do not allow loose particles to become trapped
 - a. between the concrete cylinder and the neoprene caps
 - b. or between the bearing surfaces of the extrusion controllers and the bearing blocks of the test machine.
3. Use the same surface of each neoprene cap to contact the concrete cylinders for all tests performed with that cap.
Each neoprene cap must not be used to test more than 100 cylinders.
The life of alternate pads must be verified by the agency or purchaser.

Compressive Testing Machine (cp-ct)

T 22-92, C 39-86

Purpose To identify calibration mode and type

Inspection Equipment Required

1. Calibrate load cells with meter readable to + 1%
2. Calipers, Straight Edge, Feeler Gauge

Tolerance

1. Tolerance of + 1%

Procedure

1. Apply a compressive axial load to molded cylinders or cores at a rate which is within a prescribed range until failure occurs.

As per AASHTO T 67-85 Section: C. Verification by Elastic Calibration Device.
Subsection: 13, 14 and 16-20

Scope

1. **Elastic Calibration Device**
A device for verifying the load readings of a testing machine which has:
 - (a) elastic member(s) to which loads may be applied, combined with
 - (b) mechanism or device for indicating the magnitude (or a quantity proportional to the magnitude) of deformation under load.

Concrete Test Molds (plastic) (cp-m)

ASTM C-470-94/AASHTO M 205

Purpose Inspect and test concrete molds for compliance to industry standards.

Inspection Equipment Required

See ASTM C-470.

Tolerance Molds must meet test requirements in ASTM C-470.

Inspecting new molds

1. Selection: At least three molds must be tested per shipment selected randomly.
2. Test each mold selected following ASTM C-470 procedures. If the molds tested from a particular box meet the ASTM C-470 requirements, it is assumed that all the other molds in the box are also good, and no further testing is necessary for that box.
3. If at least one of the selected molds in a box fails to meet ASTM C-470 requirements, it is assumed that several other molds in that box fail. No further testing is done, and the entire box is rejected and returned to the supplier.
4. If three boxes in the entire shipment are rejected, the whole shipment is rejected and returned to the supplier. No further testing for the shipment is necessary.

Testing New Concrete Molds

1. Molds: Meet AASHTO M-205, Section 6.
 - < Constructed in the form of right circular cylinders which stand with the cylindrical axis vertical and the top open to receive the concrete.
 - < They must be made of materials that do not react with concrete containing portland or other hydraulic cements.
 - < They must be watertight and sufficiently strong and tough to permit their use without tearing, crushing, or deforming.
2. Test Procedure:
 - a. Fill mold with water to 90 or 95 percent full.
 - b. Subject mold to tapping & jarring.
 - c. Let stand 1 hour.
 - d. Examine the mold for leaks (especially the air hole at the bottom of the mold.)
3. Examine outside the mold for curvature of the bottom. Any noticeable concave (inward curvature) or convex (outward curvature) is cause for rejection.

Reused Concrete Molds

The Materials Division does not currently reuse molds because of a high potential for problems such as water leakage, cracks, stretched diameter from extruding the first concrete cylinder under pressure, and other damage.

Conicle Molds and Tamper (cp-cm)

T84/C128

Purpose To measure critical dimensions of the sand cone and tamper used in AASHTO T84 or ASTM C128.

Inspection Equipment Required

1. Calipers or ruler readable to 1 millimeter.
2. Balance, 500 gram capacity, readable to 0.1 gram.
3. Straight edge or ruler.

Tolerance The equipment must meet the dimensions and weights as specified in the test method.

Procedure **Cone**

1. Use the calipers to measure the inside diameter at the top of the cone to the nearest 1 millimeter, taking two readings 90E apart.
Record the results.
2. Invert the cone. Repeat step #1 using the ruler.
3. Place the cone on a flat glass surface. Measure and record the depth of the cone by using the calipers and a straight edge or ruler.
Record the results.
4. Measure the thickness of the cone to the nearest 1 millimeter by taking two readings 90E apart at the top of the cone and two readings 90E apart at the bottom of the cone.
Record the results.

Tamper

1. Measure and record the diameter of the tamping face to the nearest millimeter by taking two readings 90E apart using the calipers.
Record the results.
2. Determine and record the mass of the tamper to the nearest 0.1 gram.

Direct Tension Machine (cp-dtt)

AASHTO T-314

Purpose To (a) calibration and verification of the Direct Tension Machine

Inspection Equipment Required

1. Verification standard spring for load cell and displacement transducer
2. ASTM 62 C calibrated thermometer with 0.1 EC increments

Tolerance

The Load cell, displacement transducer and temperature shall meet the accuracy requirements specified in the applicable test methods listed above.

Procedure

1. A specially designed standard (spring) shall be used to verify the operation of both the load cell and the displacement transducer. The verification standard shall have a spring rate of approximately 135 N/m and shall be capable of withstanding a 500 Newton load without damage. The manufacturer, Instron, shall supply verification standard. Instron must provide a certificate which guarantees the NIST traceability of the standard to be used. Verify calibration of the load cell and displacement transducer at least every 6 months and when measurements are suspect.
2. Verification of elongation rate- shall be verified using the verification standard (spring). The measurement shall be made at -18 EC. Allow ten minutes for the verification standard to achieve thermal equilibrium prior to running the test. To verify the elongation rate, elongation shall be plotted as a function of elapsed time. The resulting plot shall be a straight line with a slope 1.00 mm/min.
3. Temperature Detector- the temperature detector shall be calibrated as a minimum once each year. Calibration must be done at each temperature used. Direct contact between the RTD and temperature calibration device. If not agreeing, within 0.1 EC, correction applied or further calibration or maintenance must be performed. Calibrate using the same thermometer as above. Adjust the offset on the bath to correct temperature discrepancies.

4. Using the Instron “Test Specification and Procedure for doing the Verification”

located in the equipment folder, verify the following acceptance criteria:

CAcceptance Criteria for Customer's Machine

The Mean Load and Strain (maximum for each test) must fall within the following tolerances. The test is performed at 3%/min under strain control and is commanded to stop at 1% strain.

Coefficient of variation for the load must be 1.0% or less. This is as much a factor of careful specimen set-up as of actual machine variation. For each of

the following parameters, determine the difference between the mean values that were previously determined for the specimen and the values that resulted in testing the machine:

Error= Mean value from Specimen Certificate- Mean Value for machine

| | |
|----------|---|
| | Error Tolerance |
| Load | ± 1.0 N |
| Strain 1 | ± 0.0001 % (based on Gage Length of 33.8) |
| | |

All other values: Modulus, Stress, time of test, Extension, etc. should be within 6%.

- Interpretation of results

1. The load reading is outside the ± 1.0 N required by the certificate:

This may be due to a shift in calibration of the load cell, or the LVDT caused by loads induced by the method of attachment, or cantilever loads on the load cell. In this case the LVDT and Load cell must be calibrated independently without removing them from the machine. See separate test specifications for this.

2. The CV for the test is larger than 1.0%

This may be due to:

- a. The grip is not well aligned with respect to the base- use alignment disk screws to adjust.
- b. Slack in the load train- check belt tension, check that the ball screw nut is tight, check for proper preload of the ball bearing nut.
- c. Bent Verification Specimen- check for flatness

3. Time of test, or Extension is out by more than 6%
 - a. Some of the same problems listed above may apply
 - b. Check if the speed is set to 3%/min, and is under strain 1 control, set to stop at 1%
 - c. Check if the proper loadframe has been selected in Merlin setup(5525 or 55BTI)
 - d. If a 55BTI is being tested, the extension may vary by 10% from what is shown on the original test results on the certificate. Some but other values should still fall with 6%.

Ductility Apparatus (cp-da)

AASHTO T51, ASTM D113

Purpose To verify the accuracy of equipment used to perform the ductility test.

Inspection Equipment Required

1. Ruler, readable to 1 mm.
2. Timer, readable to 0.1 sec.
3. Thermometer, calibrated and readable to 0.5EC.
4. Calipers readable to 0.0001"

Tolerance

The equipment must meet the tolerances specified in test methods T 51 and D 113.

Procedure

1. For each travel speed, begin with head at zero mark and note the time to travel each 10 mm increment for a total of 150 mm.
2. Calculate travel speed for each 10 mm increment. Determine average and compare with total travel speed for 150 mm. If the calculated average and total are within 0.50 mm per minute, record the total length and time data as travel speed to the nearest 0.1 mm/min.
3. Measure with a scale the depth and height of water above and below the mold.
4. Visually inspect machine for wear on any parts, loose bolts, leaks, etc....
5. Observe and record the temperature of the water in the bath
6. Measure and record the depth of water in the bath

Molds

7. Assemble a mold.
8. Measure and record the ID of the mold.
9. With the mold assembled, measure and record:
 - a distance between centers. (Fig 1, T51)
 - b total length of the mold
 - c distance between clips
 - d shoulder
 - e width at minimum cross-section
 - f width at the mouth of the clip
 - g thickness of the mold

Dynamic Shear Rheometer (DSR)

AASHTO TP5 (DSR)

Purpose To measure the critical temperatures and torque motor response of the instrument used in AASHTO TP5.

Inspection Equipment Required

1. 2 mm thick by 25 mm diameter silicone thermistor calibrated to the nearest 0.1 EC.
2. Cannon Reference Fluid

Tolerance Equipment must meet the temperature and accuracy tolerances specified in the applicable test method.

Procedure

1. Compare temperature measurements obtained from the calibrated silicone wafer, thermistor and the DSR RTD.
Using the temperature measured by the thermistor as a reference temperature, apply an appropriate temperature correction to the temperature measurement indicated by the DSR RTD if they do not agree within $\pm 0.1\text{EC}$.
2. Measure and record the modulus and viscosity of the Cannon Reference Fluid.

Flash Cleveland Open-Cup (cp-fc)

AASHTO T48

Purpose To verify dimensional tolerances of the Flash Cleveland Open Cup.

Inspection Equipment Required

1. Calipers, readable to 0.01 mm

Tolerance The dimensions must meet the tolerances in test method T48.

Procedure

With the calipers:

1. Measure and record outside diameter of the cup at its base.
2. Measure and record the inside diameter of the cup at its base.
3. Measure and record the thickness of the bottom of the cup.
4. Measure and record from the top of the cup to the fill mark inside the cup.
5. Measure and record the outer diameter of the shoulder.

Flash Point Tag Open-Cup (cp-fp)

AASHTO T79

Purpose To verify physical dimensional tolerances of the Flash Point Tap Open Cup.

Inspection Equipment Required

1. Calipers, readable to 0.01 mm
2. Balance, 2000 g capacity, readable to 0.01 g.

Tolerance The dimensions must meet the tolerances in the test method listed above.

Procedure

With calipers:

1. Measure and record the height of the cup.
2. Measure and record the distance from the top of the cup to just under the rim on top.
3. Measure and record the diameter of the cup just under the ridge.
4. Measure and record the diameter of the cup just before it starts curving at the bottom.
5. Measure and record the depth of the cup in the center.

With balance:

6. Weigh and record the weight of the cup.

Flow Gauges (cp-fg)

AASHTO T170, T240

Purpose To verify the accuracy of the flow gauges used to control of the flow rate of air or carbon dioxide gas.

Inspection Equipment Required

1. Large volumetric flask, or other container, 1000-ml capacity or larger.
2. Gas flow meter, capable of indicating a gas flow of up to 1000 ml per min for CO₂ and 4000 ml for air.
3. Timing device readable to 0.1 s and accurate to within 0.1 s for a 60.0 s interval
4. Water bath, of sufficient size to immerse the volumetric flask.
5. Thermometer, ASTM 17C or 17F.
6. Convenient source of air or carbon dioxide gas.

Tolerance

Flow gauges must be capable of maintaining a flow rate of approximately 900 ml per min for CO₂ or 4000 ± 200 ml per min for air.

Procedure

1. Accurately measure the total volume of water in the flask.
2. Fill the both the water bath and the flask with water at 25 C ± 1E.
3. Cover the flask, invert it, and immerse it into the bath. Do not allow any air to remain in the flask.
4. Start the flow of air or carbon dioxide and set the gauge to the desired flow rate.
5. Immerse a tube connected from the flow gauge into the bath.
Start the timer at the same instant the tube outlet is place into the flask.
6. Stop the timer when all the water has been displaced from the flask.
7. Calculate the flow rate.
8. Re-set the flow gauge to a slightly higher pressure and repeat steps 2 through 7.
9. Re-set the flow gauge to a slightly lower pressure and repeat steps 2 through 7.
10. On graph paper, plot the points of actual flow versus the gauge reading from each test.
From the plot, determine the reading that indicates the desired flow.

Kinematic and Absolute Viscosity Bath (cp- kvis)

T201, T202

Purpose To verify that the operating temperature of the bath varies no more than 0.06 EC.

Inspection Equipment Required

1. NIST traceable thermometer

Tolerance Tolerances are found in test methods T201 and T202.

Procedure

1. Turn on the viscometer bath and allow it to come to test temperature
2. Leave the bath on a minimum of four hours
3. Record temperatures at half-hour intervals

LA Abrasion Machine (cp-1a)

AASHTO T96/C131

Purpose To (a) measure the critical dimensions and general operating conditions of the LA Abrasion Machine, and
(b) measure the mass of the spheres used as test charges.

Inspection Equipment Required

1. Steel rule readable to 1/16 inch (1 millimeter).
2. Stop watch readable to 0.1 second.
3. Balance with a 5 kilogram capacity readable to 1 kilogram.

Tolerance The LA machine must meet the dimensional tolerances specified in the applicable test method listed above and must be in good operating condition.
The steel spheres used to charge the LA machine must meet the mass tolerances specified in the applicable test method listed above.

Procedure **LA Machine**

1. Measure and record the inside diameter of the drum at the left and right edges to the nearest millimeter.
2. Measure and record the width and height of the opening to the nearest millimeter.
3. Measure and record the wall thickness at the left and right edge to the nearest 1/16 inch.
4. Determine if the cylinder is horizontal using a steel ball to check left to right roll.
5. Measure and record the shelf width inside the drum to the nearest 1/16 inch.
6. Measure and record the distance from the shelf to the opening in the direction of rotation.
7. Using a stop watch, determine the rpm to the nearest whole number over a five minute period; record the rpm.
8. Check that the number of revolutions is 500 by looking at the counter on the machine.

Steel Spheres

1. Determine and record the mass of each individual sphere to the nearest 1 gram.
2. Determine and record the mass of the collective charges to the nearest 1 gram.

Liquid Limit Device and Grooving Tool (cp-IIg)

T89/D4318

Purpose To inspect the physical condition and measure the critical dimensions of the liquid limit device and grooving tool.

Inspection Equipment Required

Calipers capable of measuring the critical dimensions readable to 0.001 inch.

Tolerance The critical dimensions of the liquid limit device must meet the dimensional tolerances specified in the test method.

Procedure

1. Measure and record the dimensions of the liquid limit device as indicated on figure 1 of the above test specification.
2. Inspect the liquid limit to determine that
 - a. the device is in good working order
 - b. the pin connecting the cup is not worn sufficiently to permit side play
 - c. the screws connecting the cup to the hanger arm are tight
 - d. points of contact on the cup and base are not excessively worn
 - e. the lip of the cup is not excessively worn
 - f. a groove has not been worn in the cup through long usage.

Note 1: Wear is considered excessive when

- 1) the point of contact on the cup or base exceeds approximately 13 mm (0.5 in.) in diameter, or
- 2) when any point on the rim of the cup is worn to approximately ½ the original thickness.

Although a slight groove in the center of the cup is noticeable, it is not objectionable. **If the groove becomes pronounced before other signs of wear appear, the cup should be considered excessively worn, and should be replaced.**

3. The rate of shocks of mechanically operated liquid limit devices should be similar to manually operated devices at about 2 revolutions per second.
4. With calipers, measure and record the dimensions of grooving tool and end gauge.

Manual Compaction Hammer (cp-mh)

T99/D698, T180/D1557

Purpose To measure the critical dimensions of the Proctor Hammers.

Inspection Equipment Required

1. Calipers readable to 0.001 inch.
2. Tape measure readable to 1/16 inch.
3. Balance, capacity 5 kilograms, readable to 1 gram.

Tolerance Equipment must meet the dimensional tolerances specified in the applicable test method.

Procedure

1. Using the calipers, measure the diameter of the rammer face determined by taking 2 readings 90E apart. Record the results.
2. Pull up the handle, and use the tape measure to measure the drop height of the hammer. Determine this height inside the guide sleeve. Record the results.
3. Remove the hammer from the guide sleeve; determine its mass with the balance and record the finding to the nearest 1 gram.
4. Using the calipers measure and record the diameters of the vent holes near the end of the hammer.

Marshall Hammer (cp-marh)

T-245

Purpose: To measure and compare mechanical and manual compactors.

Equipment Used:

1. Balance readable to 0.1 g
2. Tape measure of 18 inch straight edge readable to 1/16 inch
3. Calipers readable to .0001

Tolerance: As per AASHTO T-245, The difference between the two hammers should not exceed .5 lbs.

Procedure:

1. Dismantle the sliding weight. Weigh and record the weight.
2. Re-assemble the sliding weight and measure the distance from the bottom of the handle to the top of the weight. Record this distance.
3. With calipers, measure and record the diameter of the compaction face.
4. Correlate the manual and mechanical hammers every 12 months. The difference between the two hammers should not exceed 0.5 pounds.

Marshall Mechanical Compactor (cp-mm)

Purpose

Inspection Equipment Required

1. Hand operated hammer
2. Mechanical operated hammer

Tolerance Equipment must meet the dimensional tolerances specified in the applicable test method.

Procedure

1. Use the same aggregate, asphalt percentage, temperatures, and the amount of blows for steps two and three.
2. Compact marshall specimens with mechanical compactor and record the data.
3. Compact marshall specimens with a hand compactor and record the data.

Marshall Stability & Flow Test Apparatus (cp-mt)

T-245

Purpose To measure the rate of compression and the accuracy of the Marshall Stability instrument graph.

Inspection Equipment Required

1. Ruler
2. Stopwatch.
3. A load cell capable of at least 5000 lbf capacity and sensitive to 10 lbf up to 1,000 lbf and 25 lbf between 1,000 and 5,000 lbf and equipped with a micrometer dial. This compares the readings of the load cell to the rate of movement shown on the instrument graph.

Tolerance AASHTO T-245: Loading Jack with uniform vertical movement of 50.8 mm/min. and the flowmeter with gauge of .25mm divisions.

Procedure

1. Use the ruler to measure the vertical movement of loading jack and horizontal movement on the instrument graph within 30 seconds after starting the timer.

Use the stopwatch to time the amount of movement of the test head on the loading jack. It must move 2 inches in 60 seconds.

2. Use the load cell to measure the pounds of force applied.

Compare readings of the load cell to the rate of movement shown on the instrument graph to verify the accuracy of pounds reading on the graph.

Marshall Molds and Breaking Head (cp-mmbh)

T 245

Purpose To verify the dimensions of the Marshall Molds and breaking head.

Inspection Equipment Required

1. 6" Calipers capable of measuring an inside diameter of 4" and readable to 0.0001"
2. Template with a radius of 2" inner curvature

Tolerance The equipment will meet the tolerance levels specified in AASHTO T-245.

Procedure:

1. Measure and record the inside diameter of the mold to the nearest 0.0001".
2. Rotate the mold 90 degrees. Measure and record the results again.
3. Place the template into each half of the breaking head.
4. Observe any gaps, record findings.

Mechanical Compactor (cp-mc)

T99/D698, T180/D1557

Purpose To measure the critical dimensions of mechanical compactors.

Inspection Equipment Required

1. Calipers readable to 0.001 inch.
2. Tape measure readable to 1/16 inch.
3. Balance, Capacity 5 kilogram, readable to 1 gram.

Tolerance Equipment must meet the dimensional tolerances specified in the applicable test method.

Procedure

1. If circular face is used, measure with calipers and record the diameter of the rammer face by taking 2 readings 90E apart. If a sector face is used it must have an area equal to that of a circular rammer face.
2. Using the ruler, measure the drop height of hammer. Record the finding.
3. Using the balance, measure the hammer mass and record the measurement to the nearest 1 gram.

Mechanical Shaker (cp-ms)

AASHTO T11/C117, T27/C136, T88

Purpose To measure sieving's thoroughness.

Inspection Equipment Required

1. Soil or aggregate sample.
2. Balance accurate to the weights specified and of adequate capacity for the respective test methods.

Tolerance

Verify the sieving times so that the percent passing any sieve during one minute of sieving is in accordance with the value specified in the respective test procedures.

Procedure

1. Prepare soil or aggregate sample to be sieved in accordance with the respective above referenced specifications.
2. Perform sieving procedure with the mechanical shaker as many times as is customarily used for the respective test procedure.
3. Weigh and record individual sieves to determine weight of sample retained in individual sieves.
4. Use mechanical shaker to shake each individual sieve for one minute.
5. Weigh and record the individual sieves to determine the percentage of sample lost during the shaking.

Moist Cabinets, Rooms, and Water Storage Tanks (cp-mcr)

C-511-93

Purpose To verify that requirements for moist cabinets, moist rooms, and water storage tanks used for storage of paste, mortar, and concrete test specimens are being met.

Inspection Equipment Required

1. Equipment used in ASTM C-511-93.

Moist Cabinets - General

- < Maintain the atmosphere in a moist cabinet at a temperature of $73.4 \pm 3^{\circ}\text{F}$ and a relative humidity of at least 95% at all times.
- < Provide wet and dry recording thermometers for all storage units.
- < Ensure that the shelves which hold the fresh specimens are level.
- < Provide cabinets that are constructed of durable materials that have tightly fitting doors.
- < Maintain the specified relative humidity with one or more fog sprays, water sprays, or curtains of water on the inner walls that are directed so that the water will collect in a pool at or near the bottom of the moist storage section.
- < Provide automatic control of air temperature where a cabinet is located in a non air-conditioned work room, or where maintaining temperatures within the specified range is difficult.

Moist Cabinets - Monitoring Temperature and Humidity

- < Take periodic readings from the wet and dry thermometers each day with a minimum of once a day.
- < Record each reading and the time the reading is taken.
- < Determine that the readings are within the specified temperature and humidity range.
- < If the readings fall outside the acceptable range, adjust the mechanisms on the cabinet to bring the readings into the appropriate range.
- < Avoid taking humidity readings when saturation is typically below optimum, such as directly after specimens are placed in or removed from the storage unit.

Moist Room - General

- < Maintain the ambient air in a moist room at a minimum temperature of 68°F and a maximum temperature of 82.5°F , and a humidity of at least 50% at all times.

Molds - Soil Testing

T99/D698, T180/D1557, T193/D1883

Purpose To measure the critical dimensions of 4, 6, and 8 inch molds used in soil testing.

Inspection Equipment Required

1. Calipers capable of measuring the height and inside diameter of the molds readable to 0.001 inch.

Tolerance

The height and diameter of the molds measured must meet the dimensional tolerance as specified in the applicable test method listed above.

Procedure

With the calipers:

1. Measure the inside diameter of the mold to the nearest 0.001 inch. Record the results.
Rotate the mold 90E (1/4 turn) and measure the inside diameter again. Record the results.
2. Turn the mold over and repeat step one.
3. Measure the height of the mold to the nearest 0.001 inch. Record the results.
Rotate the mold 90E (1/4 turn) and measure the height again. Record the results.

Oven, General Purpose (cp-gpoven)

Purpose To verify the accuracy of dial settings on general-purpose drying ovens.

Inspection Equipment Required

1. A calibrated thermometer graduated in 1.0E C increments having a range which includes the temperature range to be checked.
2. A brass thermometer well to retain heat while the oven door is open. This is essential for a constant temperature reading.
3. A clothes pin to hold the thermometer in such a manner as to enable the operator to read the scale easily from outside or inside the oven.

Tolerance Drying ovens must be capable of maintaining a constant temperature range listed in the appropriate test methods.

Procedure

1. Place the thermometer inside the brass well with the clothes pin attached to the thermometer. Position the thermometer on the shelf where the samples are normally dried.
2. Close the oven and let it remain undisturbed. Take the first reading at least 1 hour after closing the oven.
3. Take as many readings as necessary to determine if the temperature range is within the specified tolerance. Three consecutive readings, taken no less than ½ hour apart, within tolerance allowed are adequate.
4. Adjust the temperature of the oven if an observed temperature reading is outside the tolerance specified. Allow at least ½ hour for the temperature to stabilize between each adjustment. Return to step 3.

Oven, Rolling Thin Film (cp-rtfo)

T240, D2872

Purpose To measure the critical components of a rolling thin-film oven.

Inspection Equipment Required

1. Ruler, readable to 1/8".
2. Timer, readable to 0.1 sec.

Tolerance The equipment must meet the tolerances specified in test methods T240 and D2872.

Procedure

1. Open the oven door. Measure and record the distance from the thermometer to the right side of the oven.
Measure and record the distance from the thermometer bulb to the horizontal axis running through the center of the carriage.
2. Perform the test according to T240 or D2872. After placing the test samples in the oven, record the time needed for the oven to recover to test temperature.
3. Record the number of carriage rotations in one minute.
4. Record the oven temperature at half-hour intervals until the test is complete.

Penetrometer (cp-pa)

AASHTO T49, ASTM D5

Purpose To measure, inspect, and verify the accuracy of the equipment used to perform the penetration test.

Inspection Equipment Required

1. Balance, readable to 0.01 g.
2. Microscope or eyepiece, 10X
3. Metal Block, 10.0 mm high. Metal block, 25.4 mm high.
4. Support block, 3- 3.5" high.
5. Ruler, readable to 1 mm.
6. Calibrated Stopwatch, readable to 0.1 sec.

Tolerance The equipment must meet the tolerances specified in test methods T49 and D5.

Procedure

1. Remove the spindle, 50 g and 100 g weights from the penetrometer.
Weigh and record the weight of each to the nearest 0.01 g.
2. Weigh each needle to the nearest 0.01 g.
3. Visually examine each needle with a microscope or eyepiece.
Each needle should be straight and free of burrs.
The base of each needle should be flat.
4. **If an automatic timing mechanism is used** on the penetrometer, start the calibrated stopwatch when the plunger is released and stop the calibrated stopwatch when the plunger stops.
Record the time indicated on the calibrated stopwatch to the nearest 0.1 sec.

If a manual device is used to release the plunger, check the accuracy of the timing device used over a 60 sec interval.
Record the elapsed time to the nearest 0.1sec.
5. Place the support block on the base of the penetrometer.
Place the 10 mm block on the support block.
Adjust the needle height so that its tip just touches the top of the 10 mm block.
Remove the 10 mm block and release the needle to the support block.
Adjust the instrument to measure the distance moved.

Repeat step 4 using the 25.4 mm block.
Determine dial accuracy by comparing readings with the height of the blocks.
6. Measure and record the distance from the perforated shelf to the bottom of the water bath. Measure and record the distance from the perforated shelf to the surface of the water. Measure and record the distance the thermometer is immersed in the water.
7. Observe and record the temperature of the water in the bath to the nearest 0.05 C.

Pycnometers (cp-py)

T-228

Purpose To measure the volume and the physical dimensions of the pycnometer to verify they meet tolerances.

Inspection Equipment Required

1. 7" Calipers readable to 0.0001"
2. Balance capable of weighing 2000 g. readable to 0.01 g.

Tolerance Tolerances can be found in the test methods listed above.

- Procedure**
1. Use calipers to measure and record the diameter of the stopper.
 2. Use calipers to measure and record the diameter of the hole in the stopper.
 3. Observe if the top of the stopper is smooth and plane. Record findings.
 4. Observe if the lower surface of the stopper is concave? Record findings.
 5. Measure and record the height of the concave area.
 6. Is the pycnometer capable of holding 24 to 30 mls? Record weight.
 7. Weigh the pycnometer, record weight.
 8. Calibrate the volume following the procedures in AASHTO T-228, section 6. Record results.

Rice Flasks (cp-rf)

T209

Purpose:

Inspection Equipment Required

1. Vacuum gage
2. Glass Disk
3. Astm #17F Thermometer
4. 12" Stir rod
5. Balance, capacity 5000g, readable to 0.1g

Procedure:

1. Record the dry weight of each clean, dry flask.
2. Record the weight of the glass plate to a tenth of a gram
3. Fill a 3000 ml beaker with distilled water and adjust the temperature to within 77 ± 1 E F.
4. Fill a flask with water to the top of the bottle, use the stir rod to dislodge any air bubbles trapped in the flask.
5. Hold the glass plate partially over the mouth of the flask. Pour water into the flask while sliding the glass plate in place over the mouth of the flask. There should not be any air bubbles trapped under the glass plate.
6. Dry the whole assembly begin careful to prevent air from entering the flask.
7. Record the bottle number, weight of the flask, water, and glass plate.
8. Do steps 4 through step 7 for the remainder for the flasks.

Scales and Balances (cp-sb)

M231

Purpose To verify the accuracy of scales and balances

Inspection Equipment Required

Field standard weights that comply with NIST handbook 105-1(class F)

Procedure

1. Clean any foreign material (dust, aggregate, etc.), from the scale or balance. Level the scale or balance if necessary, and adjust to zero. For electronic scales and balances, a 20 minute warm-up period is required.
2. Apply to the center of scale or balance a load equal to 1/4 capacity of the scale or balance, or as near to 1/4 capacity as weights will allow. Record the value displayed. Remove the load and record the zero-load balance.
3. Apply to the center of scale or balance, a load equal to to ½ capacity, or as near to ½ capacity as weights will allow. Record the value displayed. Remove the load and record the zero-load balance.
4. Perform the shift test at ½ capacity. Place the weights on each quarter of scale or balance (front, back, and sides singly). Record the value displayed. Remove the load and record the zero-load balance.
5. Apply to the center of the scale or balance a load equal to 3/4 capacity , or as near to 3/4 capacity as weights will allow. Record the value displayed. Remove the load and the record zero-load balance.
6. Apply to the center of scale or balance a load equal to maximum capacity, or as near to maximum capacity as weights will allow. Record the value displayed.
7. Perform the decreasing-load test. Remove the load to ½ capacity. Record the value displayed. Remove the load and record the zero-balance.
8. Compare measured values to acceptable tolerances.
9. Remove from service any scales and balances which are not in compliance, and provide for necessary repairs or replacement.

Sieves

AASHTO M92, ASTM E11

Purpose To provide instructions for inspecting and measuring the physical condition of laboratory tests sieves ranging in size from 75 millimeters (3 inches) to 0.075 millimeters (No. 200).

Inspection Equipment Required

1. A caliper readable to 0.01 millimeter (used for No. 4 and coarser)

Tolerance

Sieves must meet the physical requirements specified in AASHTO M92 (ASTM E11)

Procedure

(Steps 1 and 2 apply only to sieves having openings greater than 4.75 millimeters)

1. Select an adequate number of individual sieve openings (3 or 4) along a 45E line. Measure and record the sieve openings to verify that the size opening indicated on the label is correct.
2. Repeat step 1, rotating the sieve 90E.
3. Inspect the general condition of the sieve. Check the frame and soldered joints for cracks or holes (check for pin holes in the finer sieves).
4. Make sure the sieve has an appropriate label.
5. Check for tightness of the wires on each individual sieve

Sieve - 325

ASTM 430

Purpose To determine the correction factor for the sieve used in determining the fineness of hydraulic cement or pozzolan.

Tolerance Sieves must meet the requirements as stated in ASTM 430

Equipment required

SRM sample

#20 sieve

#325 sieve

Analytical balance capable of weighing accurately of 0.0005 grams

water spray nozzle conforming to ASTM C430 section 3.2

pressure gauge conforming to ASTM C430 section 3.3

deionized water

desiccator

Procedure

1. Open an SRM sample
Throw away 1/3 of the top sample
2. Sieve the remaining sample through a # 20 sieve
3. Isolate one sample weighing one gram. 1.0000. Put sample on the 325 sieve
4. Adjust the water pressure to read 10 ± 0.5 PSI
5. Wet sample on the sieve with approximately 50 ml of deionized water
To wet sample, begin the water spray and place the sample under it with the bottom of the spray nozzle below the top of the sieve frame (cup) about 0.5 inches. Rotate cup one revolution every second for 60 seconds.
7. Remove sample from the stream of water, and wash the sample with approximately 50 ml. Of deionized water. With a damp cloth, gently blot the underside of the screen to bring through the excess amount of water.
8. Place sieve with sample in an oven or on a hot plate heated to a temperature of 180 degrees Fahrenheit. Place the sample so there is ample room below it for the movement of air under the sieve. Too high a temperature will melt the solder holding the sieve in place and the sample will be lost.
9. Leave the sample in the oven to dry for a minimum of 45 minutes at 180E Fahrenheit or on the hot plate until the sample can move freely on the screen in all directions.
10. At the end of the drying, put the sample into a desiccator for one hour.
Weigh the sample to the nearest 0.0005 grams.
11. Calculate sieve correction factor.

$$\text{Correction factor} = \frac{(\text{reported residue (g)} - \text{found residue (g)})}{\text{Found residue (g)}} \times 100$$

Correction factor is $10 \pm$ from zero. If not in the range, recalibrate again and if it still falls outside the acceptable range, discard the sieve and replace it with a new one. Start the calibration process again with the new sieve.

12. Recalibrate the sieve after every fifth sample.

Cleaning the sieves.

Method 1:

Equipment required:

- < Ultra sonic watch bath large enough to hold the container. Low powered 150 watts maximum power input
 - < Appropriate laboratory cleaning solution. Restricted to soap or detergent-type solutions.
Do not use diluted hydrochloric acid or acetic acid to clean or rinse the sieves.
9. Place the sieve in the ultrasonic watch bath.
 10. Turn on the bath and leave it on for a complete cleaning cycle which is up to 15 minutes long.
 11. Remove the sieve and rinse with water to remove any remaining particles of lodged material.

Method 2:**Equipment needed**

- < 325 sieve
 - < Laboratory cleaning solution
 - < Hot plate
 - < Beaker or other container large enough to hold the sieve cup and cleaning solution
1. Heat the laboratory cleaning solution in the beaker over a hot plate to just below boiling.
 2. Cover the container with a watch glass to reduce evaporation, and leave the screen in the solution until the particles are loosened from the screen.
 3. Remove the screen and wash it with water.

Method 3:**Equipment needed**

- < 325 sieve
 - < Laboratory cleaning solution
 - < Beaker or other container large enough to hold the sieve cup and cleaning solution
1. Put the screen in a bath of cleaning solution and allow to remain in the solution overnight or for a period of at least 12 hours.
 2. Rinse the sieve with water, and check to see that all particles of lodged material are removed in the rinse.
 3. If particles remain in the sieve, clean the sieve again using a different cleaning method.

Slump Cones (cp-sc)

T119, C143

Purpose To measure the critical dimensions of slump cones.

Inspection Equipment Required

1. Caliper "A" capable of measuring the inside diameters of the mold. (Readable to 0.1 inch.)
2. Caliper "B" capable of measuring the thickness of the mold wall. (Readable to 0.001 inch.)
3. Ruler readable to 1/16 inch.

Tolerance The tolerance must conform to figure 1 of C-143.

Procedure

1. With calipers, measure the inside diameters of the top (smaller) end twice at right angles to each other, and average the two readings.
Record the results as the top diameter.
Follow the same procedure for the base (larger) end, and record the results.
2. Place the slump cone base on a rigid flat surface and place a ruler down the axial center of the cone so that the ruler rests on the rigid surface. Place a straight edge horizontally across the top and measure the overall height of the slump cone.
3. At three random locations at the top of the cone, and at three random locations near the bottom of the cone, measure the thickness of the metal using calipers. Average of the six readings for the mean thickness.

Sodium Sulfate Containers, Sulfate Oven (cp-ssc)

T104/C88

Purpose To visually inspect sodium sulfate containers which are used to immerse the samples of aggregate in solution, and to measure the rate of evaporation for the sulfate oven.

Inspection Equipment Required

1 liter griffin low form beakers.

Tolerance Sulfate oven should be capable of being continuously heated, and the rate of evaporation should be to the tolerances specified in the test method.

Procedure **Sodium Sulfate Containers**

1. Visually inspect the sodium sulfate containers to determine if they will permit free access of the solution to the aggregate sample, and drainage of the solution from the sample without loss of aggregate.

Sulfate Oven

1. Determine the evaporation rate from the oven by measuring the loss of water from 1 liter griffin low form beaker placed at each corner and the center of each shelf of the oven. Each beaker should initially contain 500 grams of water, and the oven will be at a temperature of $70 \pm 3^{\circ}\text{F}$ throughout the test. The test should last five hours and beakers should be weighed before and after the test.
2. The evaporation requirement applies to all test locations when the oven is empty except for the beakers of water.

Soil Hydrometers (cp-sh)

AASHTO T88/D422

Purpose Verify critical dimensions of hydrometers.

Inspection Equipment Required

1. Calipers capable of measuring the diameter of the hydrometers readable to 0.001 inches.
2. Caliper or ruler capable of measuring length of hydrometer and scale length readable to .001 inches.

Tolerance

The critical dimensions of the hydrometers should meet the dimensional tolerances specified in ASTM E100.

Procedure

1. Measure and record overall length of hydrometer.
2. Measure and record the length of the hydrometer from the tip of the bulb to the 1.0 specific gravity marking or the 0 grams per liter marking.
3. Measure and record the length of the hydrometer scale.
4. Visually inspect the outer surface of the stem and body of the hydrometer to determine that they are symmetrical about a vertical axis. There should be no abrupt changes or constrictions that would hinder thorough cleaning or tend to trap air bubbles when the instrument is immersed.
5. Visual hydrometer should always float with its axis vertical.
6. The glass should be smooth and transparent and free of bubbles, dirt or other imperfections.
7. Ballast materials should be secure in the lower part of the body with no loose material of any sort inside the hydrometer.
8. The scale should be anchored to prevent it from moving. It should show no evidence of scorching.
9. Scale should be straight without twist.

Straight Edge

T99/D698, T180/D1557

Purpose To verify the compliance of straight edge.

Inspection Equipment Required

1. Certified flat surface.
2. 0.001" Feeler Gauges

Tolerance The equipment must meet the tolerances specified in test methods referenced above.

Procedure

1. By visual inspection, determine that the straight edge is made of hardened steel, is at least 10" in length, has one beveled edge, and has at least one longitudinal surface. It should be plane within 0.1% of length within the portion used for file trimming.
2. Measure and record any space between the straight edge and the flat surface. The straight edge should not be so flexible that trimming the soil surface with the cutting edge will cause a concave soil surface.

Thermometers, Liquid-in-Glass (cp-th)

Purpose To verify the accuracy of purchased general purpose thermometers in both Fahrenheit and Centigrade measurements.

Inspection Equipment Required

1. Circulating water bath
2. NBS certified thermometer

Tolerance

General purpose thermometer must maintain an accuracy of $\pm 0.5^\circ\text{C}$. (1°F .) unless otherwise stated.

Procedure

1. Bring the water bath to a constant temperature of 25°C . (77°F .)
2. Place thermometer in the bath.
 - A. For total immersion, ensure that total immersion types are on a rack at a minimum clearance of 25 mm between the bottom of the rack and thermometers.
 - B. For partial immersion, the depth of the partial immersion should be within 1 mm ($1/32$ inch) of the immersion line on the thermometer.
3. Allow a minimum of 5 minutes for thermometer to stabilize.
4. Establish the water bath temperature.
5. Remove one thermometer at a time and read to the nearest 0.1°C (1°F).
6. Etch the date inspected on all thermometers that fall within the $\pm 0.5^\circ\text{C}$. (1°F .) tolerance.
7. Return all thermometers which are not in compliance to the supplier for replacements.

Timers (cp-ti)

T72, T201, T202, ASTM D88, D2170, D2171

Purpose To measure and verify the accuracy of timing devices.

Inspection Equipment Required

1. Atomic Clock in Boulder, Colorado
2. Timer, readable to 0.1 second, having a verified accuracy with the tolerance listed in the above test methods.

Tolerance

Timers must meet the accuracy requirements specified in the applicable test methods listed above.

Procedure

- A. Note requirement and verify
- B. Call the Atomic Clock in Boulder, Colorado (303-499-7111). At the tone start timer and note the minute given verbally. Hang up the phone.
- C. Let the timer run for at least 15 minutes.
- D. Call the Atomic Clock again and stop the timer precisely at the minute tone. Again note the minute given.
- E. Record time in minute elapsed between tones and the time elapsed on timer.

$$\begin{aligned}\% \text{ accuracy} &= (A-B)/B * 100 \\ A &= \text{reading on lab timer (sec)} \\ B &= \text{time between tones (sec)}\end{aligned}$$

Unit Weight Container - Concrete (cp-uwct)

T121, C138

Type of Equipment Unit Weight Container 1/4 cubic foot

Purpose To measure the volume of unit weight containers.

Inspection of Equipment

1. Thermometer - Capable of reading in 1E increments centigrade over a range of 0 - 50.
2. Scale or balance capable of measuring to the nearest .1 ls with an accuracy of + 1%.
3. Clear plastic plate 1/2" x 14" x 14" and plane to within .003/6".

Procedure

1. Bring all water and molds to a constant temperature of 22.7 C +/- 2 C.
2. Measure the weight or tare of the unit weight container and the plastic plate.
3. Fill the container with water at 22.7 C so that a dome is formed above the opening.
4. Place one edge of the plastic plate on one side of the unit weight container and lay the plate down on the unit weight pushing the air out to one side.
5. Ensure all the air has been evacuated and weigh the assembly. Use the following formula to calculate the volume of the vessel: $B = C - A$

- | | |
|---|--------------------------------|
| A | Tare of bucket and plate |
| B | Weight of water only |
| C | Total weight filled with water |

$$\text{Unit Weight} = B(\text{lbs.}) / 62.4 (\text{lbs./cf})$$

Unit Weight Measures - Aggregate (cp-uwa)

T19/C29

Purpose To calibrate unit weight measures.

Inspection Equipment Required

1. Plate glass 1/4 inch thick and at least 1 inch larger than the diameter of the measure.
2. Supply of water pump or chassis grease that can be placed on the rim of the container to prevent leakage.
3. Balance conforming to the requirements of AASHTO M231 for the class of general purpose.

Tolerance The equipment must meet the tolerance as specified in test method.

Procedure

1. Fill the container with water at room temperature so that a dome is formed above the opening. Place one edge of the plastic plate on one side of the measure and lay the plate down on the unit weight pushing the air out to one side.
2. Ensure that all the air is eliminated from the measure and determine the weight of water in the measure using the balance.
3. Measure the temperature of the water to determine the density and interpolate from table 3 in the specification.
4. Calculate the volume of the measure by dividing the weight of water required to fill the measure by its density.

Vacuum System (cp-vac)

T209

Purpose To measure vacuum pressure.

Inspection Equipment Required

1. Absolute pressure gauge.
2. Water vapor trap.
3. Hoses, connectors, tools, misc.

Tolerance Equipment must be capable of applying the vacuum specified in the applicable test method.
(T209)

Procedure

1. Turn on the gauge 15 minutes before beginning the test.
2. Connect the gauge to the system with the trap in line between the system and the gauge, making sure all connections are air tight.
3. Open the number of lines normally used in testing.
4. Allow the readout of the gauge to stabilize.
5. Read and record the pressure indicated on the gauge.

Water Bath (cp-wb)

T-245

Purpose To verify that the operating temperature of the bath is accurate and consistent.

Inspection equipment Required

1. Calibrated ASTM 47F or 47 C thermometer, Serial #1137

Tolerance Tolerances can be found in the test methods in T245.

Procedure

1. Turn on the water bath and allow it to come to test temperature
2. Leave it on for a minimum of four hours
3. Record temperatures at half-hour intervals

Water Bath - Penetrometer (cp-wb)

T49

Purpose To verify that the operating temperature of the bath varies no more than 0.1EC.

Inspection equipment Required

1. NIST traceable thermometer

Tolerance Tolerances can be found in the test methods in T49.

Procedure

1. Turn on the water bath and allow it to come to test temperature
2. Leave it on for a minimum of four hours
3. Record temperatures at half-hour intervals

Weighing Devices Used in the Testing of Materials

AASHTO M 231-90

| | |
|------------------|--|
| Purpose | To verify the accuracy of weighing devices. |
| Tolerance | Determined by type of balance. |
| Procedure | As procedure directs for type of balance being calibrated. AASHTO M 231-90 |
| Scope | This specification covers requirements for scales and balances in the testing of construction materials, and the loose weights used with balances. |

Calibration Record Forms

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Air Meter (cp-am)

Inspector: _____

Serial Number: _____

Location of Equipment: Central Materials Lab

Date of Calibration: _____

Calibration Interval: 12 months

Next Calibration Date: _____

Test Procedure: T-152, C-231

Calibration procedure: UDOT cp-am

Calibration Equipment Used: Threaded straight tubing, threaded bent tubing, and syringe.

Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

SITE

MEASUREMENTS

A. Temperature Verification

Temperature setting/time: Actual reading:

| | |
|-------|-------|
| _____ | _____ |
| _____ | _____ |
| _____ | _____ |

GENERAL PROCEDURE

1. Fill the base with water.
2. Screw the short piece of straight tubing into the threaded petcock hole on the underside of the cover.
Clamp cover on the base with the tube extending down into the water.
3. With both petcocks open, add water with a syringe through the petcock having the pipe extension below, until all air is forced out through the opposite petcock. Leave both petcocks open.
4. Pump up air pressure to a little beyond the predetermined initial pressure line. Wait a few seconds for compressed air to cool to normal temperature and then stabilize the gauge hand at the proper initial pressure line by pumping or bleeding off air as needed.

5. Close both petcocks and immediately press down on the thumb lever exhausting air into the base.
Wait a few seconds until the hand is stabilized.
If all the air was eliminated and the initial pressure line was correctly selected, the gauge should read 0%.
If two or more tests show a consistent variation from 0% in the result, then change initial pressure line to compensate for the variation. Use the newly established "initial pressure" line for subsequent tests.
6. Screw curved tube into the outer end of petcock.
Pressing on thumb lever and controlling flow with petcock lever, fill the 5% calibrating vessel (345 ml) level full of water from the base.
7. Release the air at the free petcock. Open the other petcock and let the water in the curved pipe run back into the base. There should be 5% air in the base.
8. With petcocks open, pump air pressure in exact manner as outlined in paragraph 4. Close petcocks and immediately press the thumb lever. Wait a few seconds for exhaust air to warm to normal temperature, and for the needle to stabilize. The dial should now read 5%.
9. If two or more consistent tests show that the gauge reads more or less than 5% air by .2% (or whatever is considered satisfactory), then remove gauge glass and reset the dial hand to 5% by turning the recalibrating screw located just below and to the right of the center diameter.
10. When gauge hand reads correctly at 5%, additional water may be withdrawn in same manner to check results at 10%, 15%, 20% ect.

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Brass Rings and Assembly (cp-br) **Date of Calibration:** _____
Inspector: _____ **Calibration Interval:** 12 months
Serial Number: _____ **Next Calibration Date:** _____
Location of Equipment: Asphalt Laboratory **Test Procedure:** T53
Calibration procedure: UDOT cp-br
Calibration Equipment and serial number: Calipers readable to 0.01 mm
Balance 2000 g capacity readable to 0.01 g.

Action Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

Check Critical Dimensions

Shouldered Rings, Fig. 1 (A):

A. O.D. at top of shoulder 22.7- 23.3 mm _____
B. I.D. at top of shoulder 19.5- 20.1 mm _____
D. I.D. at bottom of ring 15.6- 16.2 mm _____
G. Thickness of ring 6.0- 6.8 mm _____

Ring Holder, Fig. 1 (B):

B. O.D. of ring area 23.4- 24.4 mm _____
C. Total length of ring holder 75.7- 76.7 mm _____
D. I.D. of small circles 5.1- 6.1 mm _____

Ball Centering Guides, Fig. 1 (C):

D. O.D. of guide 24.3- 24.9 mm _____

Steel Balls:

Diameter 9.5 mm _____
Weight 3.45- 3.55 g _____

Comments: _____

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: CBR (cp-cbr) Molds, Annular and Slotted Weights, Penetration Piston

Inspector: _____ **Date calibrated:** _____
Serial Number: _____ **Calibration Interval:** 12 months
Location of Equipment: Soils Laboratory **Next Calibration Date:** _____
Test Procedure: T-193
Calibration procedure: UDOT cp-cbr

Calibration Equipment Used: _____

Action Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

Mold, Spacer, Tripod, Stem and Plate

Check and record dimensions shown on Figure 1.

Surcharge Weights:

Annular: Center Hole (54 mm) _____

Out Side Dia (149.2 mm) _____

Weight (2.27±0.04kg) _____

Slotted: Out Side Dia (149.2 mm) _____

Weight (2.27±0.04kg) _____

Penetration Piston:

Diameter (49.63±0.13mm) _____

Length (102+mm) _____

Comments: _____

Utah Department of Transportation
Materials Laboratory - Calibration Report

Apparatus: Compression or Loading Device (cp-ct)

Date Calibrated: _____ **Inspector:** _____

Calibration Interval: _____

Serial Number: _____ **Next Calibration Date:** _____

Location of Equipment: _____ **Calibration procedure:** UDOT cp-ct

Test Procedure: T193, T208, T216, T234, T236

Calibration Equipment Used: _____

Action Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

Temperature: _____

Scale Readings:

| | |
|---|-------|
| 1 | _____ |
| 2 | _____ |
| 3 | _____ |
| 4 | _____ |
| 5 | _____ |
| 6 | _____ |

| Recorded Load | Instrument Reading | Machine error | | Recorded Load | Instrument Reading | Machine error | |
|---------------|--------------------|---------------|---------|---------------|--------------------|---------------|---------|
| | | Lbs. | Percent | | | Lbs. | Percent |
| | | | | | | | |
| | | | | | | | |
| | | | | | | | |
| | | | | | | | |
| | | | | | | | |

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Compressive Test Machine (cp-ct) **Date of Calibration:** _____
Inspector: _____ **Calibration Interval:** 12 months

Serial Number: _____ **Next Calibration Date:** _____
Location of Equipment: Central Materials Lab **Test Procedure:** T-22-92, C-39-86
Calibration procedure: UDOT cp-ct
Calibration Equipment Used: Load Cells, Calipers, Straight Edge, Feeler Guage _____

Recommended: Repair _____ Replace _____ None _____ Other _____

| | | | |
|---------------------------------------|--------------------|------------------|---------------------|
| Bearing Block Tolerance Check: | Top | Bottom | Bottom |
| | _____ | _____ | _____ |
| | Planeness < 0.001" | Thickness > 0.9" | Other Dimensions ok |

Remarks: _____

SITE

MEASUREMENTS

A. Temperature Verification

| | |
|---------------------------|-----------------|
| Temperature setting/time: | Actual reading: |
| _____ | _____ |
| _____ | _____ |
| _____ | _____ |

General Procedure

1. Apply a compressive axial load to molded cylinders or cores at a rate which is within a prescribed range until failure occurs.

According to AASHTO T 67-85 Section: C. Verification by Elastic Calibration Device.
Subsection: 13, 14 and 16-20.

Scope:

2. Elastic Calibration Device:
A device for verifying the load readings of a testing machine which has:
 - (a) elastic member(s) to which loads may be applied, combined with
 - (b) mechanism or device for indicating the magnitude (or a quantity proportional to the magnitude) of deformation under load.

Comments: _____

Utah Department of Transportation

4501 South 2700 West, Salt Lake City, UT 84119

Materials Central Lab

Test Methods & Specifications

ASTM C-172, C-143, C-138,

C-173 or C-231, C-31, C-39

Specimen Size _____

Required 28 Day Strength _____

Concrete Compressive Strength Report

To _____ Project Engineer

Project Name & No. _____

Aggregate Source, Sand _____ Gravel _____

Air Entr. Agent _____

Concrete Class _____ Cement Source _____ Type Cement _____ Admixtures _____

Concrete Source _____ Bag Mix _____

| Technician Name | Sample Number | Name of Client & Item Cast | Cast Date | Break Date | Age | Slump | % Air | Conc. Temp | Cross Sec. Dia | Total Load | Type of Fracture | PSI | Ave PSI |
|-----------------|---------------|----------------------------|-----------|------------|-----|-------|-------|------------|----------------|------------|------------------|-----|---------|
| | | | | | | | | | | | | | |
| | | | | | | | | | | | | | |
| | | | | | | | | | | | | | |
| | | | | | | | | | | | | | |
| | | | | | | | | | | | | | |
| | | | | | | | | | | | | | |
| | | | | | | | | | | | | | |
| | | | | | | | | | | | | | |
| | | | | | | | | | | | | | |
| | | | | | | | | | | | | | |

| Org. | Acct. | Task. | CID |
|------|-------|-------|-----|
| | | | |

REMARKS: _____

Materials Engineer _____

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Concrete Test Mold (cp-mcc)

Inspector: _____

Lot Identification Number: _____

Location of Equipment: Certification Lab

Date of Calibration: _____

Calibration Interval: each box before use

Manufacturer: _____

Test Procedure: ASTM C 470-94

Calibration procedure: UDOT CP-MCC

Calibration Equipment Used: _____

Recommended: Repair _____ Replace _____ Use _____ Other _____

Remarks: Plastic re-usable mold

Certification by lot from manufacturer:

1. Plastic mold material has max water absorption of 0.5% in 24 hr, per Test Method D570.
2. Plastic mold has Izod impact toughness of 117 J/m for a 3.2 mm thick specimen, per Test Method D256.
3. After being held at -12 deg C for 24 h, the plastic must not fracture when subject to tapping and jarring per Practice C31 and C192.

Lot Test - 3 randoms samples per shipment

Limits

Pass Fail

Bottom fillet flush with sidewall?

Vert dist

3 mm

(Tolerance = ± 1.6 mm)

Horiz dist

5 mm

Resistance to damage under use (Note 1)

Summary (Note 2)

Notes:

1. Dry rodded coarse aggregate test 7.3.2, grading requirements C33, and penetration -compression testing C192.
2. Molds may be rejected for failure to comply with any provision above.

General Procedure

Molds must be constructed in the form of right circular cylinders which stand with the cylindrical axis vertical and the top open to receive the concrete. They should be made of materials that do not react with concrete containing portland or other hydraulic cements. They should be watertight and sufficiently strong and tough to permit their use without tearing, crushing, or deforming. AASHTO M-205-90 Section 6.

Fill mold with water to 90 or 95% full. Subject mold to tapping & jarring. Let stand 1 hour. Examine the mold for leaks (especially bottom hole for introducing air to remove concrete).

Examine the outside of the mold for curvature of the bottom. Any noticeable concave (inward curvature) or convex (outward curvature) is cause for rejection.

Comments: _____

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Conical Molds and Tamper (cp-cm) Date of Calibration: _____

Inspector: _____

Calibration Interval: 12 Months

Serial Number: _____

Next Calibration Date:_____

Location of Equipment: Aggregate Lab

Test Procedure: T-84

Calibration procedure: UDOT cp-cm

Calibration Equipment Used:_____

Action Recommended: Repair_____ Replace_____ None_____ Other_____

Remarks:_____

| | |
|---|--|
| Cone: | |
| Inside Dia. Top of cone (40±3 mm) 0E 90E | |
| Inside Dia. Bottom of cone (90±3 mm) 0E 90E | |
| Depth of Cone (75±3 mm) | |
| Thickness of top of Cone (0.8 mm) 0E 90E | |
| Thickness of bottom of Cone (0.8 mm) 0E 90E | |
| Tamper: | |
| Diameter of Tamper Face (25±3 mm) 0E 90E | |
| Weight of Tamper (340±15 g) | |

Comments: _____

Utah Department of Transportation

Materials Laboratory- Calibration Report

Apparatus: Direct Tension Machine Model () Date of Calibration: _____

Inspector: _____ Calibration Interval: 6 months

Serial Number: _____ Next Calibration Date: _____

Location of equipment: Asphalt Binder Laboratory Test Procedure: AASHTO T314

Calibration Procedure: cp-ddt

Calibration Equipment Used: Instron verification spring #, Calibrated ASTM

Thermometer 62 C SN

Action Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

Load cell calibrated in accordance with AASHTO T314: YES / NO

Verification of the elongation rate performed in conjunction
with the load cell and displacement transducer verification : YES / NO

(1) Measurements made at -18/C? YES / NO

(2) Elongation plotted as a function of elapsed time ? YES / NO

(3) Resulting plot a straight line with a slope of 1.00 mm/min.? YES / NO

| Instron Verification Specimen Serial # | | Mean/CV | Verification Results | Difference | Tolerance |
|---|--|---------|-------------------------|------------|-----------|
| | | | | | |
| Load at Max Strain 1, N | | Mean | | | +/- 1.0 N |
| | | CV | | | < 1.0% |
| Modulus, Mpa | | Mean | | | < 6 % |
| CV of the 20 sec loop | | CV | | | < 6 % |

Verification of the temperature detector performed daily by comparing the RTD output with a calibrated liquid-in-glass thermometer

(1) Calibrated at each temperature used? YES / NO

(2) Direct contact between RTD and temperature calibration device? YES / NO

(3) If not agreeing within " 0.1/C, correction applied or further
calibration or maintenance performed? YES / NO

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Dynamic Shear Rheometer (SHRP) (cp-dsr) Date of Calibration: _____

Inspector: _____

Calibration Interval: 6 Months

Serial Number: _____

Next Calibration Date: _____

Location of Equipment: Asphalt Laboratory

Test Procedure: T51

Calibration Procedure: cp-dsr

Calibration Equipment _____

Action Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

1. Specimen Trimmer:

A. Straight Edge Width: Minimum 4 mm _____ mm.

B. Physical Inspection: Any rough edges ? Yes / No

11. Base Test Plates:

25 mm Plate:

A. Diameter 25 mm plate 25.00 ± 0.05 mm

| | | | |
|-------------|-------------|-------------|-------------------|
| 1) _____ mm | 2) _____ mm | 3) _____ mm | Average: _____ mm |
|-------------|-------------|-------------|-------------------|

B. Height of Raised Portion 2.00 to 5.00 mm

| | | | |
|-------------|-------------|-------------|-------------------|
| 1) _____ mm | 2) _____ mm | 3) _____ mm | Average: _____ mm |
|-------------|-------------|-------------|-------------------|

C. Physical Inspection

Scratches or Imperfections? Yes / No

8 mm plate:

A. Diameter 8 mm plate 8.00 ± 0.05 mm

| | | | |
|-------------|-------------|-------------|-------------------|
| 1) _____ mm | 2) _____ mm | 3) _____ mm | Average: _____ mm |
|-------------|-------------|-------------|-------------------|

B. Height of Raised Portion 2.00 to 5.00 mm

| | | | |
|-------------|-------------|-------------|-------------------|
| 1) _____ mm | 2) _____ mm | 3) _____ mm | Average: _____ mm |
|-------------|-------------|-------------|-------------------|

C. Physical Inspection

Scratches or Imperfections? Yes/ No

111. Spindles:

25 mm Spindle:

A. Diameter 25 mm plate 25.00 ± 0.05 mm

| | | | |
|-------------|-------------|-------------|-------------------|
| 1) _____ mm | 2) _____ mm | 3) _____ mm | Average: _____ mm |
|-------------|-------------|-------------|-------------------|

B. Physical Inspection

Scratches or Imperfections? Yes/ No

8 mm Spindle:

A. Diameter 8 mm plate 8.00 ± 0.05 mm

| | | | |
|-------|-------|-------|-------------|
| 1) mm | 2) mm | 3) mm | Average: mm |
|-------|-------|-------|-------------|

B. Physical Inspection

Scratches or Imperfections? Yes/ No

Gap:

A. Measure the gap between the base plate and spindle.

8 mm gap is.....2 mm ± 0.001 mm Yes/ No

25mm gap is.....1 mm ± 0.001 mm Yes/ No

1V. Thermistor: S/N_____

A. Thickness # 2.00 mm _____ mm.

B. Thermistor calibration (See attached worksheet)

C. Run Temperature Calibration for Software based on
Calibrated Thermistor Yes/ No

D. Physical Inspection
Wafer surface not bent and wires OK ? Yes/ No

V. Verification of Overall Calibration:

A. Bohlin representative's yearly DSR/Torque check.
(See attached calibration report)

B. Utilizing a Cannon Viscosity Standard # _____ at;

| Temperature, | EC Cannon Viscosity Standard mPa s (cP) | DSR Calculated Viscosity mPa s (cP) | DSR Tolerance $\pm 1\%$ of Cannon Standard |
|--------------|--|---|--|
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |

Note: Conversion from Pascal to Centipoise:

$(G^*(\text{Pascal})/\text{Oscillating Frequency (radians/second)}) \times 1000$

Comments: _____

1V. Thermistor Calibration continued:

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Ductility Apparatus (cp-da), molds Date of Calibration: _____

Inspector: _____ **Calibration Interval:** 12 months

Serial Number: 20457-16 **Next Calibration Date:**

Location of Equipment: Asphalt Laboratory **Test Procedure:** T51

Calibration procedure: UDOT cp-da

Calibration Equipment Used: Calipers readable to 0.0001"

Action Recommended: Repair_____ Replace_____ None_____ Other_____

Remarks:

| SITE | SPECIFIED | MEASURED |
|--|------------------------|----------|
| A. Travel Speed | 47.5 to 52.5 mm/min | _____ |
| B. Functions without undue Vibration | | _____ |
| C. Space for at least 25 mm water above and below mold | | _____ |
| D. Physical Inspection | | _____ |
| E. Bath Temperature | 4 plus or minus 0.05 C | _____ |
| F. Depth of water in Bath | greater than 50 mm | _____ |
| G. Thickness of mold | 9.9 to 10.1 mm | _____ |
| H. Width of mold at minimum cross section | 9.9 to 10.1 mm | _____ |

Comments:

General Procedure

1. For each travel speed, begin with head at zero mark and note the time to travel each 10 mm increment for a total of 150 mm.
2. Calculate travel speed for each 10 mm increment. Determine average and compare with total travel speed for 150 mm. If the calculated average and total are within 0.50 mm per minute, record the total length and time data as travel speed to the nearest 0.1 mm/min.
3. Measure with a scale the depth and height of water above and below the mold.
4. Visually inspect machine for wear on any parts, loose bolts, leaks, etc....

Ductility Apparatus (cp-da), molds Page 2

- Observe and record the temperature of the water in the bath
- Measure and record the depth of water in the bath

Molds

7. Assemble a mold.
8. Measure and record the ID of the mold.
9. With the mold assembled, measure and record:
 - a distance between centers. (Fig 1, T51) _____
 - b total length of the mold _____
 - c distance between clips _____
 - d shoulder _____
 - e width at minimum cross-section _____
 - f width at the mouth of the clip _____
 - g thickness of the mold _____

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Flash Cleveland Open-Cup (cp-fc)

Date of Calibration: _____

Calibration Interval: 12 Months

Next Calibration Date:_____

Test Procedure: T48

Calibration procedure: UDOT cp-fc

Calibration Equipment Used: CALIPERS READABLE TO 0.01 mm SER #

Action Recommended: Repair_____ Replace_____ None_____ Other _____

Remarks: _____

| | # | # | # |
|---|-------|-------|-------|
| 1. O.D. at the base of cup 67.5-69.0 mm | _____ | _____ | _____ |

2. I.D. at the base of cup 62.4-64.0 mm _____

3. Thickness of bottom 2.8-3.6 mm _____

4. From the top of cup to fill mark 9-10 mm _____

5. O.D. of flange (J) 97-101 mm _____

6. Total depth of cup 32.5-34 mm _____

Comments: _____

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Flash Point Tag Open-Cup (cp-fp) **Date of Calibration:** _____
Inspector: _____ **Calibration Interval:** 12 months
Serial Number: _____ **Next Calibration Date:** _____
Location of Equipment: Asphalt Laboratory **Test Procedure:** T79
Calibration procedure: UDOT cp-fp
Calibration Equipment Used: Calipers, readable to 0.0001" ser # _____
Balance, readable to 0.01g ser# _____

Action Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

| | # 1 | # 2 | # 3 |
|--|-------|-------|-------|
| 1. Height of cup 50.0- 53.2 mm | _____ | _____ | _____ |
| 2. Rim to top of cup 7.1-8.7 mm | _____ | _____ | _____ |
| 3. Diameter under the ridge 53.1-56.9 mm | _____ | _____ | _____ |
| 4. Diameter before curvature 49.2-52.4 | _____ | _____ | _____ |
| 5. Inside depth of cup 47.6 mm | _____ | _____ | _____ |
| 6. Weight of cup is not more than 95 g | _____ | _____ | _____ |

Comments: _____

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Flow Gauges (cp-fg)

Date of Calibration: _____

Inspector: _____

Calibration Interval: 12 months

Serial Number: _____

Next Calibration Date: _____

Location of Equipment: _____

Test Procedure: T170, T240

Calibration procedure: UDOT cp-fp

Calibration Equipment Used: Timing device, ser # _____
ASTM 17C or 17F Thermometer, Ser # _____

Action Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

| | #1 | # | 2 |
|---|-------|-------|-------|
| | | | |
| 1. Total volume of water at 25 ± 1 E C in the filled flask. | _____ | _____ | _____ |
| 2. Time required to evacuate water. | _____ | _____ | _____ |
| 3. Flow rate = Volume m./Time, min. | _____ | _____ | _____ |
| 4. Plot actual flow versus the gauge reading. | | | |
| 5. Reading from plot for desired flow: | | | _____ |

Utah Department of Transportation

Materials Laboratory - Calibration Report

Kinematic Viscometers

Apparatus: Zeitfuchs Cross- Arm

Date of Calibration: _____

Inspector: _____

Calibration Interval: 36 months

Serial Number: _____

Next Calibration Date: _____

Location of Equipment: Asphalt Laboratory

Test Procedure: T201, A3

Use Calibration Form: UDOT cp-vb

Calibration Equipment Used:

Viscosity Standards

Calibrated Thermometer

Calibrated Timer

Action Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

MEASURED

A. Efflux time (t) _____

B. Viscosity of Standard (v) _____

Use formula $C = v/t$ _____

C. Constant C _____

Comments: _____

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: LA Abrasion Machine (cp-la) **Date of Calibration:** _____
Inspector: _____ **Calibration Interval:** 24 Months
Make: Soiltest **Model:** _____ **Serial Number:** UDOT 59-44-78
Next Calibration Date: _____ **Test Procedure:** AASHTO T-96
Location of Equipment: Aggregate Lab. **Calibration Procedure:** UDOT CP-LA
Calibration Equipment Used: _____

Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

| | |
|--|--|
| Inside Diameter(711+/-5mm) | |
| Inside Length(508+/-5mm) | |
| Opening Dimensions | |
| Axis Horizontal(within 1 in 100 slope) | |
| Removable steel shelf full length projecting inward(89+/-2mm) | |
| Shelf thickness (mm) | |
| Shelf position at least greater than 1.27m from opening on outside diameter of cylinder in direction of rotation | |
| Shelf not bent | |
| Check for backlash or slip in driving mechanism | |
| Charge (steel balls all weigh between 390 and 445g) Record weights (to nearest g) | |
| 1 | |
| 2 | |
| 3 | |
| 4 | |
| 5 | |
| 6 | |
| 7 | |

| | |
|--|--|
| 8 | |
| 9 | |
| 10 | |
| 11 | |
| Total | |
| Determine rpm of drum (30 to 33) | |
| Check that counter correctly counts to 500 | |

Comments: _____

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Liquid Limit Device and Grooving Tool (cp-llg)

Date of Calibration: _____

Inspector: _____

Calibration Interval: 12 months

Serial Number: _____

Next Calibration Date: _____

Location of Equipment: Geotechnical Lab

Test Procedure: T-89

Calibration procedure: UDOT cp-llg

Calibration Equipment Used: Calipers readable to 0.001 in.

Action Recommended: Repair_____ Replace_____ None_____ Other _____

Remarks: _____

| | |
|--|--|
| Measure dimensions of the liquid limit device and grooving tools and record on attached Figure 1 | |
| Rate of shocks (should be about 2 per second) | |
| Does pin connecting cup permit excessive side play? | |
| Screws connecting cup to hanger are tight? | |
| Point of contact on cup and base excessively worn? (see Note 1) | |
| Lip of cup excessively worn? (see Note 1) | |
| Excessive groove worn in cup? (see Note 1) | |
| Is device in good working order? | |

Note 1: Wear is considered excessive when

- 1) the point of contact on the cup or base exceeds approximately 13 mm (0.5 in.) in diameter, or
- 2) when any point on the rim of the cup is worn to approximately ½ the original thickness.

Although a slight groove in the center of the cup is noticeable, it is not objectionable.

If the groove becomes pronounced before other signs of wear appear, the cup should be considered excessively worn, and should be replaced.

Liquid Limit Device and Grooving Tool (cp-llg) Page 2

A base which is excessively worn may be refinished as long as

- 1) the thickness does not exceed the tolerance shown in Figure 1 by more than -2.5 mm (-0.1 in.)
- 2) the distance between the cup at the cam follower and the base is maintained within the

tolerances specified in Figure 1.

| | | | | |
|-----------------|-----------|--|--|--|
| Gage End | 10±.2mm | | | |
| Cutting Edge | 2.0±.1mm | | | |
| Width | 13.5±.1mm | | | |
| Depth of Groove | 10±.2mm | | | |

Comments: _____

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Manual Compaction Hammer (cp-mh) **Date of Calibration:** _____
Inspector: _____ **Calibration Interval:** 12 months
Serial Number: _____ **Next Calibration Date:** _____
Location of Equipment: Geotechnical Lab **Test Procedure:** T-99 and T-180
Calibration procedure: UDOT cp-mh

Calibration Equipment Used: _____

Action Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

| | |
|---|--|
| Rammer Fall Height (304.8 ± 1.524 mm) | |
| Rammer Fall Height (457.2 ± 1.525 mm) | |
| Rammer Circular Face Diameter (50.8 ± 0.254 mm): 0 degrees 90 degrees | |
| If not circular, area of face (2006 to 2047 mm ²) | |
| Weight of Hammer (2.495 ± 0.009 kg) | |
| Weight of Hammer (4.536 ± 0.009 kg) | |
| 4 Vent Holes (9.5+ mm) | |

Comments: _____

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Marshall Hammer (cp-marh)

Date of Calibration: _____

Inspector: _____

Calibration Interval: 12 months

Serial Number: N/A

Next Calibration Date: _____

Manufacturer: Humboldt

Model: N/A

Location of Equipment: _____

Test Procedure: AASHTO T-245

Calibration procedure: UDOT cp-marh

Calibration Equipment Used:

Balance readable to 0.1g

Tape measure or 18" straight edge readable to 1/16" Calipers
readable to .0001"

Action Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

| | | | |
|--------------|--------------|------------------|-----------|
| AASHTO SPEC. | 98.4mm | 4527g. - 4545.0g | 18 inches |
| Hammer I.D. | Dia. of face | Sliding Weight | Free fall |

Is correlation with a manual hammer within .5 lbs.? _____

Is correlation done yearly? _____

Procedure:

1. Dismantle the sliding weight. Weigh and record the weight.
2. Re-assemble the sliding weight and measure the distance from the bottom of the handle to the top of the weight. Record this distance.
3. With calipers, measure and record the diameter of the compaction face.
4. Correlation between the manual and mechanical hammers is done every 12 months. The difference between the two hammers should not exceed 0.5 lbs.

Comments: _____

Utah Department of Transportation

Materials Laboratory - Calibration Report

Marshall Mechanical Compactor (cp-mm)

Inspector: _____

Serial Number: N/A

Manufacturer: Humboldt

Location of Equipment: Bituminous Lab

Date of Calibration: _____

Calibration Interval: 36 months

Next Calibration Date: _____

Model: H-1356

Test Procedure: UDOT cp-mm

Calibration Equipment Used: Hand operated hammer, mechanical operated hammer

Action Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

Mechanical bulk density: _____ **Hand bulk density:** _____

Procedure:

1. Use the same aggregate, asphalt percentage, temperatures, and the amount of blows for steps two and three.
2. Compact marshall specimens with mechanical compactor and record the data.
3. Compact marshall specimens with a hand compactor and record the data.

Comments: _____

Utah Department of Transportation

Materials Laboratory - Calibration Report

Marshall Molds & Breaking Head (cp-mmbh) **Date of Calibration:** _____
Inspector: _____ **Calibration Interval:** 12 months
Serial Number: N/A **Next Calibration Date:** _____
Manufacturer: Humboldt **Model:** H-1342
Location of Equipment: Bituminous Lab **Test Procedure:** AASHTO T-245

Calibration Equipment Used: 1) 6" calipers capable of measuring an inside diameter of 4" and readable to 0.0001",
2) template with a radius of 2" inner curvature

Action Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

Diameter specification 3.995" to 4.005"

| Top of mold | | Bottom of mold | |
|-------------|----------|----------------|----------|
| 1. _____ | 2. _____ | 1. _____ | 2. _____ |
| 1. _____ | 2. _____ | 1. _____ | 2. _____ |
| 1. _____ | 2. _____ | 1. _____ | 2. _____ |
| 1. _____ | 2. _____ | 1. _____ | 2. _____ |
| 1. _____ | 2. _____ | 1. _____ | 2. _____ |
| 1. _____ | 2. _____ | 1. _____ | 2. _____ |
| 1. _____ | 2. _____ | 1. _____ | 2. _____ |
| 1. _____ | 2. _____ | 1. _____ | 2. _____ |
| 1. _____ | 2. _____ | 1. _____ | 2. _____ |

Any gaps between template and breaking head? _____
Does the breaking head assembly bind? _____

Procedure:

1. Measure and record the inside diameter of the mold to the nearest 0.0001".
2. Rotate the mold 90 degrees. Measure and record the results again.
2. Place the template into each half of the breaking head.
3. Observe any gaps, record findings.

Comments: _____

Utah Department of Transportation

Materials Laboratory - Calibration Report

Marshall Stability & Flow Test Apparatus (cp-mt) **Date of Calibration:** _____

Inspector: _____

Calibration Interval: 12 months

Serial Number: St#44-1037

Next Calibration Date: _____

Manufacturer: Humboldt

Model: H-1339

Location of Equipment: Bituminous Lab

Test Procedure: AASHTO T-245

Calibration procedure: UDOT cp-mt

Calibration Equipment Used: Stopwatch

Measuring tape or 12" ruler

Load cell with a sensitivity of 10 lbf up to 1000 lbf (44.5 N to 4.45 kN) & from 1000 to 5000 lbf (4.45 kN to 22.2 kN) sensitivity.

Action Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

Distance moved: _____ **Time:** 1 minute

Rate of movement: _____ (AASHTO standard - 2in/min or 50.8 mm/min)

Load at 1,000 lbf: _____ **Load at 5,000 lbf:** _____

Procedure:

1. Use the ruler to measure the vertical movement of loading jack and horizontal movement on the instrument graph within 30 seconds after starting the timer.

Use the stopwatch to time the amount of movement of the test head on the loading jack. It must move 2 inches in 60 seconds.

2. Use the load cell to measure the pounds of force applied.

Compare readings of the load cell to the rate of movement shown on the instrument graph to verify the accuracy of pounds reading on the graph.

Comments: _____

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Mechanical Compactor (cp-mc)

Date of Calibration: _____

Inspector: _____

Calibration Interval: 12 months

Serial Number: _____

Next Calibration Date: _____

Location of Equipment: _____

Test Procedure: T99/D698, T180/D1557

Calibration procedure: UDOT cp-mc

Calibration Equipment Used: _____

Action Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

Rammer Fall Height (304.8 ± 1.524 mm) _____

Rammer Fall Height (457.2 ± 1.525 mm) _____

Rammer Circular Face Diameter (50.8 ± 0.254 mm) _____

0 degrees

90 degrees

If not circular, area of face (2006 to 2047 mm²) _____

Weight of Hammer (2.495 ± 0.009 kg) _____

Weight of Hammer (4.536 ± 0.009 kg) _____

Comments: _____

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: **Mechanical Shaker (cp-ms)**

Date of Calibration: _____

Inspector: _____

Calibration Interval: _____

Serial Number: _____

Next Calibration Date: _____

Location of Equipment: _____

Test Procedure: _____

Calibration procedure: UDOT cp-ms

Calibration Equipment Used: _____

Action Recommended: Repair_____ Replace_____ None_____ Other _____

Remarks: _____

| Sieve Size | Sieve Weight (g) | Sieve+Soil Weight (g) After Mechanical Shaking | Sieve+Soil Weight (g) After Hand Shaking 1 Min. | Percent Loss During Shaking |
|------------|------------------|--|---|-----------------------------|
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |

Comments: _____

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Moist Cabinet, Room (cp-mcr)

Date of Calibration: _____

Inspector: _____

Calibration Interval: 12 months

Serial Number: _____

Next Calibration Date: _____

Location of Equipment: Central Materials Lab

Test Procedure: C-511-93

Calibration Equipment Used: Calibrated thermometer and humidity measuring device.

Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

A. Temperature Verification Temperature setting/time: Actual reading:

| | |
|-------|-------|
| _____ | _____ |
| _____ | _____ |
| _____ | _____ |

Moist Cabinets - General

- < Maintain the atmosphere in a moist cabinet at a temperature of $73.4 \pm 3\text{EF}$ and a relative humidity of at least 95% at all times.
- < Equip all storage units with wet and dry recording thermometers.
- < Shelves on which fresh specimens are placed should be level.
- < Provide cabinets which are constructed of durable materials with tightly fitting doors.
- < Maintain the specified relative humidity with one or more fog sprays, water sprays, or curtains of water on the inner walls that are directed so that the water will collect in a pool at or near the bottom of the moist storage section.
- < Provide automatic control of air temperature where a cabinet is located in a non air-conditioned work room, or where maintaining temperatures within the specified range is difficult.

Moist Cabinets - Monitoring Temperature and Humidity

- < Take periodic readings from the wet and dry thermometers each day with a minimum of once a day.
- < Record each reading and the time the reading is taken.
- < Determine that the readings are within the specified temperature and humidity range.
- < If the readings fall outside the acceptable range, adjust the mechanisms on the cabinet to bring the readings into the appropriate range.
- < Avoid taking humidity readings when saturation is typically below optimum, such as directly after specimens are placed in or removed from the storage unit.

Moist Room - General

- < Maintain the ambient air in a moist room at a minimum temperature of 68E F and a maximum temperature of 82.5EF , and a humidity of at least 50% at all times.

Comments: _____

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Mold -soils testing, compaction (cp-m) **Date of Calibration:** _____ **Inspector:** _____

Calibration Interval: 12 months

Serial Number: _____ **Next Calibration Date:** _____

Location of Equipment: Geotechnical Lab

T e s t

Procedure: T-99, T-180, T-193

Calibration procedure: UDOT cp-m

Calibration Equipment Used: _____

Action Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

| Dimension Measured | Measurement |
|---|---------------|
| Inside Diameter 4 inch mold (4+/-0.016 in) | |
| | At 90 degrees |
| Invert 4 inch mold diameter | |
| | At 90 degrees |
| Height of 4 inch mold (4.584+/-0.005 in) | |
| | At 90 degrees |
| Volume of 4 inch mold (0.0333+/-0.0003 cu ft) | |
| Inside Diameter 6 inch mold (6+/-0.026 in) | |
| | At 90 degrees |
| Invert 6 inch mold diameter | |
| | At 90 degrees |
| Height of 6 inch mold (4.584+/-0.005 in) | |
| | At 90 degrees |
| Volume of 6 inch mold (0.075+/-0.00075 cu ft) | |

Comments: _____

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Oven, General Purpose (cp-gpoven) **Date of Calibration:** _____
Inspector: _____ **Calibration Interval:** 4 Months
Serial Number: _____ **Next Calibration Date:** _____
Location of Equipment: _____ **Calibration procedure:** UDOT cp-gpoven

Calibration Equipment Used: Calibrated Thermometer, 1E increments, Ser # _____
Clothes pin
Brass thermometer well
Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

Site Measurements

A. Temperature Verification

| Temperature setting/time: | Actual reading: |
|---------------------------|-----------------|
| _____ | _____ |
| _____ | _____ |
| _____ | _____ |

General Procedure

1. Place the thermometer inside the brass well with the clothes pin attached to the thermometer. Position the thermometer on the shelf where the samples are normally placed.
2. Close the oven and let it remain undisturbed. Take the first reading at least 1 hour after closing the oven.
3. Take as many readings as necessary to determine if the temperature range is within the specified tolerance. Three consecutive readings, taken no less than ½ hour apart, within tolerance allowed, are adequate.
4. Adjust the temperature of the oven if an observed temperature reading is outside the tolerance specified. Allow at least ½ hour for the temperature to stabilize between each adjustment. Return to step three.

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Oven, Rolling Thin Film (cp-rtfo) **Date of Calibration:** _____
Inspector: _____ **Calibration Interval:** 12 months
Serial Number: _____ **Next Calibration Date:** _____
Location of Equipment: Asphalt Laboratory **Test Procedure:** T240
 Calibration procedure: UDOT cp-rtfo
Calibration Equipment Used: _____

Action Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

| SITE | SPECIFIED | MEASURED |
|------------------------------|--------------------|----------|
| A. Thermometer placement | 50.8 plus or minus | _____ |
| Vertical | 3.2 mm | _____ |
| B. Thermometer placement | 25 mm | _____ |
| C. Test temperature recovery | 10 minutes max. | _____ |
| D. Carriage rotation | 15 plus or minus | _____ |
| | 0.2 R/min | _____ |
| E. Oven temperatures | 163 plus or minus | _____ |
| | 0.5 C | _____ |

Procedure

1. Open the oven door. Measure and record the distance from the thermometer to the right side of the oven. Measure and record the distance from the thermometer bulb to the horizontal axis running through the center of the carriage.
2. Perform the test according to T240 or D2872. After placing the test samples in the oven, record the time needed for the oven to recover to test temperature.
3. Record the number of carriage rotations in one minute.
4. Record the oven temperature at half-hour intervals until the test is complete.

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Penetrometer (cp-pa) **Date of Calibration:** _____
Inspector: _____ **Calibration Interval:** 6 months
Serial Number: _____ **Next Calibration Date:** _____
Location of Equipment: Asphalt Laboratory **Test Procedure:** T49
 Calibration procedure: UDOT cp-pa
Calibration Equipment Used: _____

Action Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

| SITE | SPECIFIED | MEASURED |
|--|----------------------------|----------|
| A. 50 g weight | 50.0 plus or minus 0.05 g | _____ |
| 100 g weight | 100.0 plus or minus 0.05 g | _____ |
| needle and spindle | 47.5 plus or minus 0.05 g | _____ |
| B. Needle weight | _____ | |
| C. Physical inspection of needle | _____ | |
| D. Timing mechanism | _____ | |
| E. Dial accuracy | _____ | |
| F. Dist. perforated shelf to bottom of water bath greater than 50 mm | | _____ |
| Dist. perforated shelf to water surface greater than 100 mm | | _____ |
| Dist. thermometer immersed in water 150 plus or minus 15 mm | | _____ |
| G. Water bath temperature | | _____ |

Comments: _____

Page 3-99

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Rice Flasks (cp-rf)
Inspector: _____
Serial Number: N/A
Manufacturer: Various
Location of Equipment: Bituminous lab

Date of Calibration: _____
Calibration Interval: 12 months
Next Calibration Date: _____
Model: _____
Test Procedure: AASHTO T-209
Calibration procedure: UDOT cp-rf

Calibration Equipment Used: Vacuum gage, Glass disk, Astm #17F Thermometer, 12" Stir rod, Balance capacity 5000g, readable to 0.1g

Action Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

20 inches of mercury minimum? _____
#200 screen on ends of stoppers? _____
Gauge placed at end of system? _____

| Flask I.D. | Dry weight of flask | Flask filled with water and glass plate |
|------------|---------------------|---|
|------------|---------------------|---|

| | | |
|-------|-------|-------|
| _____ | _____ | _____ |
| _____ | _____ | _____ |
| _____ | _____ | _____ |
| _____ | _____ | _____ |

Weight of glass plate _____

GENERAL PROCEDURE

1. Record the dry weight of each clean, dry flask.
2. Record the weight of the glass plate to a tenth of a gram
3. Fill a 3000 ml beaker with distilled water and adjust the temperature to within 77 ± 1 E F.
4. Fill a flask with water to the top of the bottle, use the stir rod to dislodge any air bubbles trapped in the flask.
5. Hold the glass plate partially over the mouth of the flask. Pour water into the flask while sliding the glass plate in place over the mouth of the flask. There should not be any air bubbles trapped under the glass plate.
6. Dry the whole assembly begin careful to prevent air from entering the flask.
7. Record the bottle number, weight of the flask, water, and glass plate.
8. Do steps 4 through step 7 for the remainder for the flasks.

Comments: _____

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Sand Equivalent Equipment
Inspector: _____

Date of Calibration: _____

Serial Number: _____

Calibration Interval: 12 months

Location of Equipment: Aggregate Lab.

Next Calibration Date: _____

Calibration Equipment Used: _____

Test Procedure: T-176 (Standard)

Action Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

| | |
|---|--|
| Graduated Plastic Cylinder, rubber stopper, irrigator weighted foot assembly and syphon assembly | Measure and record dimensions shown on attached Figure 1, show measurements on Figure 1. |
| 1 gallon bottle of working calcium chloride solution shelf height above working surface? (36 in. +/- 1 in.) | |
| 3 oz. tin box: Diameter (2.25 in.) Capacity (85 +/- 5 ml) | |
| 4 in. wide mouth funnel | |
| Timer: a b $(a-b)/b \times 100$ | |
| Mechanical Shaker (Soil Test-232, No. 144): Throw 8 ± 0.4 in. Cycles per minute (175 ± 2) | |
| Mechanical Kelvinator Thermo-Resistive Unit: Measure mho's in degrees Kelvin | |

Comments: _____

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Saybolt Viscometer
Inspector: _____
Serial Number: _____
Location of Equipment: Asphalt Laboratory

Date of Calibration: _____
Calibration Interval: 36 Months
Next Calibration Date: _____
Test Procedure: T72, section 9
Use Calibration form for UDOT cp-vb

Calibration Equipment Used: Calibrated Thermometer
Viscosity Standard
Calibrated Timer

Action Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

| SITE | SPECIFIED | MEASURED OR CERTIFIED |
|-------------------------------|-----------|-----------------------|
| A. Laboratory Temperature | 20-30E C | _____ |
| B. Viscosity Bath Temperature | 50.0E C | _____ |
| C. Viscosity of Standard (V) | > 90 sec | _____ |
| D. Efflux time (t) | | _____ |

E. Does efflux time of Viscosity Standard differ from Certified

Saybolt viscosity value by more than 0.2 percent? yes _____ no _____

F. If yes, Calculate a correction factor , (F) using $F=V/t$ _____

Note: Do not use viscometers or orifices in referee testing if they require corrections greater than 1.0 percent.

Comments: _____

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Scales and Balances (cp-sb) Date of Calibration: _____
Inspector: _____ Calibration Interval: _____
Serial Number: _____ Next Calibration Date: _____
Manufacturer: _____ Model: _____
Location of Equipment: _____ Calibration Procedure: UDOT cp-sb
Standard/Weights Identification: _____

Action Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

SCALE AND BALANCE WORKSHEET

Scale Capacity: _____ Acceptance Tolerance: _____

Value of Scale Division: _____ Maintenance Tolerance: _____

| | Measured | | Zero Load | |
|-------------------------|----------|-------|-----------|---------|
| Increasing - Load Test: | Load | Value | Tolerance | Balance |
| 1/4 Capacity | _____ | _____ | _____ | _____ |
| 1/2 Capacity | _____ | _____ | _____ | _____ |
| Shift Test | _____ | _____ | _____ | _____ |
| 3/4 Capacity | _____ | _____ | _____ | _____ |
| Full Capacity | _____ | _____ | _____ | _____ |
| Decreasing - Load Test | _____ | _____ | _____ | _____ |
| Zero - Load Balance | _____ | _____ | _____ | _____ |

Within Tolerance: _____ Out of Tolerance: _____

Comments: _____

Utah Department of Transportation

Materials Laboratory - Calibration Report

Sample Printout of Results

Apparatus: Sieves AASHTO M-92/ASTM E-11/Utah UDOT (Std. Spec. for Wire Cloth Sieves)
Project: R4 Main-Rich 81MATRL4/SFD/4266/8732/6219
Feature: 12 new full/ Warehouse) to Larry Gay 896-13-61337KED
Manufacturer: W.S. Tyler/USA No. 99759825
Sieve Size: 2.360 (No. 8) Sieve No. 002789.Z#
By: JAL (UDOT) Date: 9-9-97

| No. | Aperture Dimensions | | Wire Diameter Readings | |
|-----|---------------------|--------------|------------------------|--------------|
| | Shoot (mm) | Warp (mm) | Shoot (mm) | Warp (mm) |
| 1 | 2.35800 | 2.35700 | 1.06700 | .99000 |
| 2 | 2.32700 | 2.31800 | 1.07100 | 1.00600 |
| 3 | 2.34300 | 2.39400 | 1.05400 | .99900 |
| 4 | 2.28500 | 2.39400 | 1.04900 | .99200 |
| 5 | 2.35900 | 2.38800 | 1.03500 | .99400 |
| 6 | 2.32500 | 2.38300 | 1.05700 | .99800 |
| 7 | 2.34900 | 2.29900 | 1.05000 | 1.00100 |
| 8 | 2.34900 | 2.41700 | 1.06400 | .99200 |
| 9 | 2.35800 | 2.38600 | 1.06300 | .99300 |
| 10 | 2.34300 | 2.34700 | 1.05500 | .99000 |
| 11 | 2.31800 | 2.39400 | 1.05100 | 1.00300 |
| 12 | 2.38500 | 2.33200 | 1.05900 | 1.00600 |
| 13 | 2.33600 | 2.39000 | 1.04800 | .99700 |
| 14 | 2.35400 | 2.38100 | 1.04600 | .98700 |
| 15 | 2.35200 | 2.40500 | 1.05600 | .99100 |
| 16 | 2.33400 | 2.34100 | 1.06300 | .99600 |
| 17 | 2.33100 | 2.41800 | 1.05600 | 1.00100 |
| 18 | 2.36500 | 2.31200 | 1.05700 | .99700 |
| 19 | 2.35900 | 2.40600 | 1.05500 | 1.00300 |
| 20 | 2.30700 | 2.26400 | 1.04500 | 1.00500 |

Statistical Calculations

| | | | | |
|---------------|---------|---------|---------|--------|
| Mean | 2.34685 | 2.36710 | 1.05495 | .99705 |
| Std Dev | .02017 | .04366 | .00827 | .00583 |
| Cond Limit | 2.28000 | 2.44000 | 1.15000 | .85000 |
| Req. #reading | 1 | 6 | 0 | 0 |

Aperture Calculations

| | | |
|---|---|---------|
| Average actual shoot aperture (limits are 2.280 to 2.440mm) | = | 2.35675 |
| Average actual warp aperture (limits are 2,280 to 2,440mm) | = | 2.36710 |
| Percent aperture greater than 2,525 mm (limit 5% max.) | = | 0.00000 |
| Largest actual aperture (limit = 2.600 mm) | = | 2.41800 |

Wire Diameter Calculations

| | | |
|---|---|---------|
| Average actual diameter shoot (limits are 0.850 to 1.150mm) | = | 1.05495 |
| Average actual diameter warp (limts are 0.850 to 1.150mm) | = | .99705 |

Visual Check OK **Specification check** OK

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Sieve 325 (cp-sv325)
Inspector: R. Zimmerman
Serial Number: _____
Manufacturer: _____
Location of Equipment: _____

Date of Calibration: _____
Calibration Interval: _____
Next Calibration Date: _____
Model: _____
Test Procedure: ASTM 430
Calibration Procedure: UDOT cp-sv325

Equipment required:

SRM sample, #20 sieve, #325 sieve.

Analytical balance capable of weighing accurately of 0.0005 grams.

Water spray nozzle conforming to ASTM C320 section 3.2.

Pressure gauge conforming to ASTM C430 section 3.3.

Dionized water, desiccator.

Sample test data from logbook pages maintained in the Cement Lab.

| <i>Cup #</i> | <i>Factor</i> |
|--|---------------|
| 1 .0830-.0808 = .0022 ($.0022 \div .080$)(100) = | 2.7228 |
| 2 .0830-.0794 = .0036 ($.0036 \div .0794$)(100)= | 4.5340 |
| 3 .0830-.0848 = -.0018 ($-.0018 \div .0848$)(100)= | -2.1226 |
| 4 .0830-.0832 = -.0002 ($-.0002 \div .0832$)(100)= | -0.2404 |
| 5 .0830-.0846 = 0.0015 ($-.0015 \div .0846$)(100)= | -1.7730 |

| Sample # | Cup # | D | C | Used | % Return | Run Date |
|--------------|-------|-------|---------|------|--------------------------|----------|
| F-1 | 1 | .1856 | 2.7228 | 1/ | 19.1 | 1.16.97 |
| F-2 | 2 | .2302 | 4.5340 | 1/ | 24.1 | 1.16.97 |
| F-3 | 3 | .1992 | -2.1226 | 1/ | 19.5 | 1.16.97 |
| F-4 | 4 | .1829 | -0.2404 | 1/ | 18.2 | 1.16.97 |
| F-5 | 5 | .1706 | -1.7730 | 1/ | 16.8 | 1.16.97 |
| Fed 123 | 1 | .0523 | 2.7228 | 2/ | 94.6278 Passing 94.63 | 1.31.97 |
| Fed 124 | 2 | .0422 | 4.5340 | 2/ | 95.5364 Passing 95.54 | 1.31.97 |
| 97-R4-M-0021 | 1 | .1983 | 2.7228 | 2/ | 20.3 | 3.12.97 |
| F-1 | 1 | .1784 | 2.7228 | 3/ | 18.3 | 5-5-97 |

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Slump Cone (cp-sc)

Date of Calibration: _____

Inspector: _____ **Calibration Interval:** 12 months
Serial Number: H-3640 **Next Calibration Date:** _____
Location of Equipment: Central Materials Lab **Test Procedure:** T-119, C-143
Calibration Procedure: UDOT cp-sc

Calibration Equipment Used: Caliper A readable to 0.1 inch (diameters)
Caliper B readable to 0.0001 inch. (wall thickness)
Ruler readable to 1/16 inch.

Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

| | Base Diameter 8 ± 1/8 " | Wall Thickness 0.045 " Minimum | |
|---------|-----------------------------------|--|---------------|
| | | Top | Bottom |
| Average | _____ | _____ | _____ |
| | _____ | _____ | _____ |
| | _____ | _____ | _____ |
| | Top Diameter 4 ± 1/8 " | Height of Mold 12 ± 1/8 " | |
| Average | _____ | _____ | |
| | _____ | | |
| | _____ | | |

Comments: _____

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Sodium Sulfate Containers (cp-ssc), Sulfate Oven

Date of Calibration: _____ **Calibration procedure:** cp-ssc

Inspector: _____ **Calibration Interval:** 12 Months

Serial Number: _____ **Next Calibration Date:** _____

Location of Equipment: Aggregate Lab. **Test Procedure:** T-104

Calibration Procedure: UDOT cp-ssc

Calibration Equipment Used: _____

Action Recommended: Repair _____ Replace _____ None _____ Other _____

| Beaker # | #1 | #2 | #3 | #4 | #5 |
|----------------------------------|-------|-------|-------|-------|-------|
| Tare & H₂O | | | | | |
| Tare | | | | | |
| Wt H₂O | 500.0 | 500.0 | 500.0 | 500.0 | 500.0 |

Time _____

| | | | | | |
|----------------------------------|--|--|--|--|--|
| Tare & H₂O | | | | | |
| Tare | | | | | |
| Wt H₂O | | | | | |

Time _____

| | | | | | |
|------------------------------|-------|-------|-------|-------|-------|
| Initial Wt. H ₂ O | 500.0 | 500.0 | 500.0 | 500.0 | 500.0 |
| Final Wt. H ₂ O | | | | | |
| Wt loss | | | | | |

| | | | | | |
|------------------------------|-------|-------|-------|-------|-------|
| Wt loss | | | | | |
| Initial Wt. H ₂ O | 500.0 | 500.0 | 500.0 | 500.0 | 500.0 |

| | | | | | |
|---------------------|--|--|--|--|--|
| Percent Loss | | | | | |
|---------------------|--|--|--|--|--|

A minimum of 20% loss is required **Oven Temp.** _____ (70±3F) _____

Sodium Sulfate Containers adequate to permit free access of the solution to the aggregate sample and drainage of the solution for the sample without loss of aggregate? _____

Comments: _____

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Soil Hydrometer (cp-sh)

Inspector: _____

Serial Number: _____

Location of Equipment: Soils Laboratory

Date of Calibration: _____

Calibration Interval: 24 months

Next Calibration Date: _____

Test Procedure: T-88

Calibration procedure: cp-sh

Calibration Equipment Used: _____

Action Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

| | |
|--|--|
| Length of Hydrometer (278 to 280 mm) | |
| Length for tip to 1.0 sp.g. or 0 g/l marking (244 to 246 mm) | |
| Length of scale (82 to 84 mm) | |
| Symmetrical about vertical axis? | |
| Abrupt changes/constructions that could trap air bubbles? | |
| Floats with axis vertical? | |
| Glass smooth, transparent and free of bubbles and imperfections? | |
| Ballast material secure? Any loose material in hydrometer? | |
| Scale anchored without twist and not scorched? | |

Comments: _____

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Straight Edge (cp-se)

Inspector:

Serial Number: _____

Location of Equipment: Soils Laboratory

Date of Calibration: _____

Calibration Interval: 6 months

Next Calibration Date: _____

Test Procedure: T99, T180

Calibration procedure: cp-se

Calibration Equipment Used: Certified flat surface
0.001" Feeler Gauges

Action Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

Straight Edge OK _____ YES _____ NO

Inspection procedure:

1. By visual inspection, determine that the straight edge is made of hardened steel, is at least 10" in length, has one beveled edge, and has at least one longitudinal surface. It should be plane within 0.1% of length within the portion used for file trimming.
2. Measure and record any space between the straight edge and the flat surface. The straight edge should not be so flexible that trimming the soil surface with the cutting edge will cause a concave soil surface.

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Thermometers (cp-th)

Date of Calibration: _____

Inspector: _____

Calibration Interval: 6 months

Serial Number: _____

Next Calibration Date: _____

Location of Equipment: _____

Test Procedure: T201,T202,T228,T49,T51

Calibration procedure: cp-th

Calibration Equipment Used

Serial Number:

Action Recommended: Repair_____ Replace_____ None_____ Other _____

Remarks: _____

| Laboratory | Therm.type,ser# | | Reading at 25EC (77FF) | | True temp. | Correction | Remarks |
|------------|-----------------|--|------------------------|--|------------|------------|---------|
| | | | | | 25.0 (77) | | |
| | | | | | 25.0 (77) | | |
| | | | | | 25.0 (77) | | |
| | | | | | 25.0 (77) | | |
| | | | | | 25.0 (77) | | |
| | | | | | 25.0 (77) | | |
| | | | | | 25.0 (77) | | |
| | | | | | 25.0 (77) | | |
| | | | | | 25.0 (77) | | |
| | | | | | 25.0 (77) | | |
| | | | | | 25.0 (77) | | |
| | | | | | 25.0 (77) | | |
| | | | | | 25.0 (77) | | |

If + add corrections to thermometer to get true temperature.

If - subtract correction from thermometer to get true temperature.

Thermometers (cp-th), page 2

| Laboratory | Therm type. ser.# | Reading @50 C | True temp. | Correction | Remarks |
|------------|-------------------|---------------|---------------|------------|---------|
| | | | 50.0 C (122F) | | |
| | | | 50.0 C (122F) | | |
| | | | 50.0 C (122F) | | |

| Laboratory | Therm type. ser.# | Reading @60 C | True temp. | Correction | Remarks |
|------------|-------------------|---------------|---------------|------------|---------|
| | | | 60.0 C (140F) | | |
| | | | 60.0 C (140F) | | |
| | | | 60.0 C (140F) | | |
| | | | 60.0 C (140F) | | |
| | | | 60.0 C (140F) | | |
| | | | 60.0 C (140F) | | |
| | | | 60.0 C (140F) | | |
| | | | 60.0 C (140F) | | |

| Laboratory | Therm type. Ser# | Reading @163 C | True temp. | Correction | Remarks |
|------------|------------------|----------------|----------------|------------|---------|
| | | | 163.0 C (325F) | | |
| | | | 163.0 C (325F) | | |
| | | | 163.0 C (325F) | | |
| | | | 163.0 C (325F) | | |
| | | | 163.0 C (325F) | | |

If + add corrections to thermometer to get true temperature.

If - subtract correction from thermometer reading to get true temperature.

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Thermometer, Liquid-in-Glass (cp-th) **Date of Calibration:** _____
Inspector: _____ **Calibration Interval:** 6 months
Serial Number: _____ **Next Calibration Date:** _____
Range: _____ F or C in _____ Division **Calibration Procedure:** UDOT cp-th
Immersion: _____ **Mark:** _____ ASTM _____
Location of Equipment: _____
Calibration Equipment Used: _____

Action Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

| Certified Temperature | Thermometer Reading | Correction |
|-----------------------|---------------------|------------|
| _____ | _____ | _____ |
| _____ | _____ | _____ |
| _____ | _____ | _____ |
| _____ | _____ | _____ |

Tolerance: General purpose thermometers must maintain an accuracy of ± 1 F unless otherwise stated.

GENERAL PROCEDURE

1. Bring the water bath to a constant temperature of 85 F.
2. Place thermometer in the bath.
 - A. For total immersion, ensure that total immersion types are on a rack at a minimum clearance of 1" between the bottom of the rack and thermometers.
 - B. For partial immersion, the depth of the partial immersion must be within 1 mm of the immersion line on the thermometer.
3. Allow a minimum of 5 minutes for thermometer to stabilize.
4. Establish the water bath temperature.
5. Remove one thermometer at a time and read to the nearest 1 F.
6. Etch the date inspected on all thermometers that fall within the ± 1 F tolerance.
7. Return all thermometers which are not in compliance to the supplier for replacements.

Reference: National Institute of Standards and Technology Test No. G33783. This thermometer has been tested by comparison with standards certified by NIST. If the correction is + the true temperature is higher than the reading. If the correction is - the true temperature is lower than the thermometer reading. All temperatures are based on IPTS-68. If the ice point is included, subsequent change in its reading will change all other readings by the same amount.

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Timers (cp-ti) **Date of Calibration:** _____
Inspector: _____ **Calibration Interval:** 6 months
Serial Number: _____ **Next Calibration Date:** _____
Location of Equipment: _____ **Calibration Procedure:** UDOT cp-ti
Calibration Equipment Used: Atomic Clock- Boulder Colorado
Action Recommended: Repair _____ Replace _____ None _____ Other _____
Remarks: _____

| SITE | SPECIFIED | MEASURED |
|-------------------------|--|----------|
| A. Graduation intervals | 0.1 second | _____ |
| B. % accuracy | Accurate to within 0.05% when tested over intervals of not less than 15 minute | _____ |

Comments: _____

GENERAL PROCEDURE (atomic clock accessible)

- Note requirement and verify
- Call the Atomic Clock in Boulder, Colorado (303-499-7111). At the tone start timer and note the minute given verbally. Hang up the phone.
- Let the timer run for at least 15 minutes.
- Call the Atomic Clock again and stop the timer precisely at the minute tone. Again note the minute given.
- Record time in minute elapsed between tones and the time elapsed on timer.

% accuracy= (A-B)/B *100
A= reading on lab timer (sec)
B= time between tones (sec)

Comments: _____

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Unit Weight Container - Concrete (cp-uwct)
of Calibration: _____

Date _____

Inspector: _____

Calibration Interval: 12 months

Serial Number: H-3661

Next Calibration Date: _____

Location of Equipment: Central Materials Lab

Test Procedure: T-121, C-138

Calibration procedure: UDOT cp-uwct

Calibration Equipment Used: Calibrated thermometer, scale or balance, and clear plastic plate.

Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

SITE

MEASUREMENTS

A. Temperature Verification

Temperature setting/time:

Actual reading:

Comments: _____

GENERAL PROCEDURE

1. Bring all water and molds to a constant temperature of $22.7^{\circ}\text{C} + 2^{\circ}\text{C}$.
2. Measure the weight or tare of the unit weight container and the plastic plate.
3. Fill the container with water at 22.7°C so that a dome is formed above the opening.
4. Place one edge of the plastic plate on one side of the unit weight container and lay the plate down on the unit weight pushing the air out to one side.
5. Ensure all the air has been evacuated and weigh the assembly. Use the following formula to calculate the volume of the vessel: $c-A + \text{Volume in cu ft } 62.4$
 - A Tare of bucket and plate
 - B Weight of water only
 - C Total weight filled with water

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Unit Weight Measures - Aggregate (cp-uwa)

Date of Calibration: _____

Inspector: _____ **Calibration Interval:** 12 months

Serial Number: _____ **Next Calibration Date:** _____

Location of Equipment: Aggregate Lab

Test Procedure: T-19

Calibration procedure: UDOT cp-uwa

Calibration Equipment Used: _____

Action Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

| | |
|--|--|
| Temperature of Water and Molds, EC | |
| Weight of Unit Weight Container and Plastic Plate - a, kg | |
| Weight of Unit Weight Container, Plastic Plate, and full of water at - c, kg | |
| Volume of Unit Weight Container, m ³ (c-a)/b b=unit weight of water kg/m ³ for given temperature from Table 3 in Procedure T-19 | |

Comments: _____

Utah Department of Transportation

Materials Laboratory - Calibration Report

Vacuum Capillary Viscometers

Inspector: _____

Serial Number: _____

Location of Equipment: Asphalt Laboratory

Date of Calibration: _____

Calibration Interval: _____

Next Calibration Date: _____

Test Procedure: T202.X4

Use Calibration Form :UDOT cp-vb

Calibration Equipment Used: Viscosity Standards
Calibrated Stopwatch
Calibrated Thermometer

Action Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

| | Timing Marks | Trial 1 | Trial 2 | |
|---|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| A. Time (t) between timing marks. | _____ _____ _____ | _____ _____ _____ | _____ _____ _____ | |
| B. Viscosity of viscosity Standard at calibration Temperature (v) | | | STANDARD _____ | |
| | | K | K | K |
| | Bulb | Trial 1 | Trial 2 | Average |
| C. Calculate viscometer bulb calibration factor (K) | _____ _____ _____ _____ | _____ _____ _____ _____ | _____ _____ _____ _____ | _____ _____ _____ _____ |

Note: the duplicate determinations of calibration constant, K, for each bulb must agree within 2 per cent of their mean. Use the formula: $K = v/t$

General Procedure: see AASHTO T 202, X4.3 through X4.3.1.8

Comments: _____

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Vacuum System (cp-vac) **Date of Calibration:** _____
Inspector: _____ **Calibration Interval:** 12 months
Serial Number: St#44-1122 **Next Calibration Date:** _____
Manufacturer: Welch-Duo-Seal **Model:** 14502 L-05
Location of Equipment: _____ **Test Procedure:** In house procedure/T209
Calibration procedure: UDOT cp-vac

Calibration Equipment Used: Vacuum gauge capable of 30 Hg

Action Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

Reading in psig: _____

Reading in Hg: _____

Procedure:

1. Turn on the gauge 15 minutes before beginning the test.
2. Connect the gauge to the system with the trap in line between the system and the gauge, making sure all connections are air tight.
3. Open the number of lines normally used in testing.
4. Allow the readout of the gauge to stabilize.
5. Read and record the pressure indicated on the gauge.

Comments: _____

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Viscometer Bath (cp-vb, cp-kvis) **Date of Calibration:** _____
Inspector: _____ **Calibration Interval:** 6 months
Serial Number: _____ **Next Calibration Date:** _____
Manufacturer: _____ **Model:** _____
Location of Equipment: Asphalt Lab **Test Procedure:** ASHTO T-245
Calibration procedures: UDOT cp-vb,
UDOT cp-kvis

Calibration Equipment Used: _____

Action Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

| Time | Temperature |
|-------|-------------|
| _____ | _____ |
| _____ | _____ |
| _____ | _____ |
| _____ | _____ |

Bath temperature maintained to within $\pm 0.06^\circ\text{C}$ of test temperature for operating purposes.

Yes _____ No _____

Comments: _____

Procedure:

1. Turn on the viscometer bath and allow it to come to test temperature. Leave on for a minimum of four hours.
2. Record temperatures at half hour intervals.

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Water Bath (cp-wb)

Inspector: _____

Serial Number: _____

Manufacturer: VWR Scientific

Location of Equipment: Bituminous Lab

Date of Calibration: _____

Calibration Interval: 6 months

Next Calibration Date: _____

Model: 1245-PC

Test Procedure: ASHTO T-245

Calibration procedure: UDOT cp-wb

Calibration Equipment Used: Calibrated ASTM 47F or 47C thermometer Serial # 1137

Action Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

| Time | Temperature |
|-------|-------------|
| _____ | _____ |
| _____ | _____ |
| _____ | _____ |
| _____ | _____ |
| _____ | _____ |
| _____ | _____ |
| _____ | _____ |

Procedure:

1. Turn on the water bath and allow it to come to test temperature. Leave on for a minimum of four hours.
2. Record temperatures at half hour intervals.

Comments: _____

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Water Bath - Penetrometer (cp-wb) **Date of Calibration:** _____
Inspector: _____ **Calibration Interval:** 6 months
Serial Number: _____ **Next Calibration Date:** _____
Manufacturer: _____ **Model:** _____
Location of Equipment: Asphalt Lab **Test Procedure:** ASHTO T-49
Calibration procedure: UDOT cp-wb

Calibration Equipment Used: _____

Action Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

| Time | Temperature |
|-------|-------------|
| _____ | _____ |
| _____ | _____ |
| _____ | _____ |
| _____ | _____ |
| _____ | _____ |
| _____ | _____ |
| _____ | _____ |

Procedure:

Bath temperature maintained to within $\pm 0.1^{\circ}\text{C}$. Of the test temperature for operating purposes

YES _____ NO _____

Comments: _____

Utah Department of Transportation

Materials Laboratory - Calibration Report

Apparatus: Weighing Devices - Materials Tests (cp-wd)

Date of Calibration: _____ Calibration Interval: _____
Inspector: _____ Next Calibration Date: _____
Serial Number: _____ Model: _____
Manufacturer: _____ Test Procedure: ASHTO M 231-90
Location of Equipment: _____ Calibration procedure: UDOT cp-wd

Calibration Equipment Used: _____

Action Recommended: Repair _____ Replace _____ None _____ Other _____

Remarks: _____

Scale And Balance Worksheet

Scale Capacity: _____ Acceptance Tolerance: _____

Value of Scale Division: _____ Maintenance Tolerance: _____

| | Measured | | Zero Load | |
|-------------------------|----------|-------|-----------|---------|
| Increasing - Load Test: | Load | Value | Tolerance | Balance |
| 1/4 Capacity | _____ | _____ | _____ | _____ |
| 1/2 Capacity | _____ | _____ | _____ | _____ |
| Shift Test | _____ | _____ | _____ | _____ |
| 3/4 Capacity | _____ | _____ | _____ | _____ |
| Full Capacity | _____ | _____ | _____ | _____ |
| Decreasing - Load Test | _____ | _____ | _____ | _____ |
| Zero - Load Balance | _____ | _____ | _____ | _____ |

Within Tolerance: _____ Out of Tolerance: _____

Comments: _____